



Примеры хроматограмм ZB-5HT Inferno

Butter Triglycerides



Column: Zebron ZB-5HT Inferno
Dimensions: 15 meter x 0.32 mm x 0.10 μ m
Part No.: 7EM-G015-02
Injection: On-Column @ 103 °C, 2 μ L
Carrier Gas: Helium @ 1.8 mL/min (constant flow)
Oven Program: 100 °C to 400 °C @ 14 °C/min for 10 min
Detector: FID @ 400 °C
Sample: Butter

Mineral Oil Using H-53 Conditions




Column: Zebron ZB-5HT Inferno
Dimensions: 30 meter x 0.25 mm x 0.10 μ m
Part No.: 7HG-G015-02
Injection: On-Column @ 53 °C, 0.1 μ L
Oven Program: 50 °C for 6 min to 400 °C @ 20 °C/min for 15 min
Carrier Gas: Helium @ 1.3 mL/min (constant flow)
Detector: FID @ 415 °C
Sample: Fuel was 10 mg/mL in dichloromethane with 50 ppm markers
1. Decane (C10)
2. Tetracontane (C40)

Polybrominated Diphenyl Ethers (PBDEs)


Column: Zebron ZB-5HT Inferno™
Dimensions: 15 meter x 0.25 mm x 0.10 µm
Part No.: 7EG-G015-02
Injection: On-Column @ 73 °C, 0.5 µL
Carrier Gas: Helium @ 1.5 mL/min (constant flow)
Oven Program: 70 °C to 160 °C at 25 °C/min to 350 °C at 10 °C/min for 5 min
Detector: ECD @ 400 °C
Sample: Sample was 2.5 ppm in Isooctane

- | | |
|-------------|-------------|
| 1. BDE-25 | 11. BDE-99 |
| 2. BDE-28 | 12. BDE-85 |
| 3. BDE-75 | 13. BDE-154 |
| 4. BDE-49 | 14. BDE-153 |
| 5. BDE-71 | 15. BDE-138 |
| 6. BDE-47 | 16. BDE-183 |
| 7. BDE-66 | 17. BDE-190 |
| 8. BDE-77 | 18. BDE-203 |
| 9. BDE-100 | 19. BDE-209 |
| 10. BDE-119 | |



Fast H53 Analysis of Diesel Fuel


Column: Zebron ZB-5HT Inferno™
Dimensions: 15 meter x 0.32 mm x 0.10 µm
Part No.: 7EM-G015-02
Injection: On-Column @ 63 °C, 0.1 µL
Carrier Gas: Helium @ 2.7 mL/min (constant flow)
Oven Program: 60 °C to 375 °C @ 25 °C/min
Detector: FID @ 400 °C
Sample: Diesel Fuel was 200 ppm in Dichloromethane with Internal Standards at 50 ppm
 1. Decane (C10)
 2. Tetracontane (C40)



**Free and Total Glycerin in B-100 Biodiesel
Methyl Esters: ASTM Method D 6584**

Column: Zebtron ZB-5HT Inferno
Dimensions: 15 meter x 0.32 mm x 0.10 µm
Part No.: 7EM-G015-02
Injection: On-Column @ 53 °C, 1 µL
Carrier Gas: Helium @ 3 mL/min (constant flow)
Oven Program: 50 °C for 1 min to 180 °C @ 15 °C/min to 230 °C
@ 7 °C/min to 380 °C @ 30 °C/min for 10 min
Detector: FID @ 380 °C
Sample: A 2.0 m x 0.53 mm guard column was connected to the analytical column per ASTM method requirements

| | |
|------------------------------|---------------------|
| 1. Glycerol | 6. Tricarpin (IS 2) |
| 2. Butanetriol (IS 1) | 7. Diglycerides |
| 3. Ester | 8. 1,3-Diolein |
| 4. Monoglycerides | 9. Triglycerides |
| 5. 1-Monooleoyl-rac-glycerol | 10. Triolein |

**Free and Total Glycerin in Biodiesel:
Method EN 14105**

Column: Zebtron ZB-5HT Inferno
Dimensions: 15 meter x 0.32 mm x 0.10 µm
Part No.: 7EM-G015-02
Injection: On-Column @ 50 °C, 1 µL
Carrier Gas: Helium @ 5.79 cm/sec (constant flow)
Oven Program: 50 °C for 1 min to 180 °C @ 15 °C/min to 230 °C
@ 7 °C/min to 370 °C @ 10 °C/min (hold 5 min)
Detector: FID @ 380 °C
Sample: A 5 m x 0.53 mm guard column was connected to the analytical column

| | |
|-----------------------|---------------------|
| 1. Glycerol | 4. Tricarpin (IS 2) |
| 2. Butanetriol (IS 1) | 5. Diolein |
| 3. Monoolein | 6. Triolein |

