

Carbohydrate and Organic Acid Analysis

- Excellent resolution and column-to-column reproducibility
- Easy, accurate quantitation from sharper peak shapes
- Longer column lifetimes and faster run time capability from lower backpressures
- Baseline separation of critical sample components due to higher efficiencies

Rezex ion-exclusion HPLC columns achieve reproducible, accurate separations based on multiple modes of interaction. Available in 4 % and 8 % cross-linked sulfonated styrene-divinylbenzene (SDVB) and multiple ionic forms (calcium, sodium, hydrogen, potassium, lead, and silver) for a wide range of selectivities. USP L17, L19, L22, L34, and L58 packings available.



If you are not completely satisfied with the performance of any Rezex column, as compared to a competing product of the same size and phase, simply return the Rezex column with your comparative data within 45 days for a FULL REFUND.

Use Rezex for carbohydrate, oligosaccharide, and organic acid separations:

- Drug formulation and excipient analysis
- Food and beverage quality control testing
- Fermentation reaction monitoring and recovery testing for biofuels


Recommended alternative to Bio-Rad® Aminex®, Supelco® SUPELCOGEL™, and Waters® Sugar-Pak™ (see p. 250)

Find the Column For Your Application

Phases Available	Description	Applications	Additional Notes
RCM-Monosaccharide (L19 packing)*	8 % cross-linked resin CALCIUM ionic form	Monosaccharides and sugar alcohols from sweeteners and corn and cane sugars Di, tri, and tetra saccharides	– Our most commonly used column type – Easy regeneration with calcium nitrate solutions
RHM-Monosaccharide (L17 packing)*	8 % cross-linked resin HYDROGEN ionic form	Monosaccharides in combination with organic acids, fatty acids, alcohols, ketones, neutral compounds, or inorganic salts	– Versatile column, generally run with a mobile phase of deionized water
RAM-Carbohydrate	8 % cross-linked resin SILVER ionic form	Selectivity complementary to other Rezex column types	
RSO-Oligosaccharide	4 % cross-linked resin SILVER ionic form	High resolution of oligosaccharides up to 18 degrees of polymerization (Dp)	– Guard column is recommended to protect the ionic integrity of the matrix
RNO-Oligosaccharide	4 % cross-linked resin SODIUM ionic form	High resolution of oligosaccharides	
RPM-Monosaccharide (L34 packing)*	8 % cross-linked resin LEAD ionic form	Monosaccharides and sugar alcohol analysis. Cellobiose, glucose, xylose, arabinose, mannose and other cellulose products	
RNM-Carbohydrate (L58 packing)*	8 % cross-linked resin SODIUM ionic form	For matrices which contain high concentration of inorganic sodium, i.e. molasses	– Easily regenerated to the original ionic strength
ROA-Organic Acid (L22 packing)*	8 % cross-linked resin HYDROGEN ionic form	Organic acids alone or in combination with carbohydrates, alcohols, fatty acids, or neutral compounds; Amino sugars; ethanol, acetic acid, glycerol, and standard alcohol mixtures	– Selectivity can be altered by changing the pH as well as the type of dilute mineral acid used as the mobile phase
RFQ-Fast Acid	8 % cross-linked resin HYDROGEN ionic form	Rapid screening of fruit quality; ethanol, acetic acid, glycerol, and standard alcohol mixtures	– Analytes are routinely chromatographed under 5 minutes
RKP-Potassium	8 % cross-linked resin POTASSIUM ionic form	Analysis of glyphosate	
RCU-USP Sugar Alcohols (L19 packing)*	8 % cross-linked resin CALCIUM ionic form	For sugar analysis according to the USP procedures	– Sorbitol and mannitol can be resolved using simple isocratic conditions

* United States Pharmacopeia (USP)






Bioethanol Fermentation Monitoring


- Easy quantitation of ethanol fermentation broth components
- Monitor starches, sugars, organic acids, and ethanol in one run
- Reliable lactic acid and acetic acid monitoring
- Increase throughput by reducing run times 50 % with 150 mm column length

Monitoring the key reaction components throughout the fermentation process is crucial for maximizing ethanol recovery. Rezex ROA is uniquely suited for the separation and analysis of simple and complex sugars, organic acids, and ethanol within a fermentation broth sample. With results easily obtained through an isocratic run, Rezex ROA is instrumental in helping you to accurately determine what critical steps need to be taken to ensure the maximum yield is achieved during your fermentation run.

Rezex ROA has the ability to achieve excellent baseline separation between Dp3 and Dp4+, which have proven to be a challenge within the bioethanol industry. It is this great baseline separation that affords scientists the opportunity to utilize a shorter column dimension. By using the 150 x 7.8 mm Rezex ROA column, you are able to decrease the run time by 50 % when compared to the average run time on a 300 x 7.8 mm column.



Column: Rezex ROA-Organic Acid
Dimensions: 300 x 7.8 mm
Part No.: 00H-0138-K0
Mobile Phase: 0.005 N Sulfuric Acid
Flow Rate: 0.6 mL/min
Detection: RI @ 40 °C
Temperature: 60 °C
System: Shimadzu® Prominence® LC-20A System
Sample: 1. Dp4+ 5. Lactic Acid
 2. Dp3 6. Glycerol
 3. Maltose 7. Acetic Acid
 4. Glucose 8. Ethanol



Column: Rezex ROA-Organic Acid
Dimensions: 150 x 7.8 mm
Part No.: 00F-0138-K0
Mobile Phase: 0.005 N Sulfuric Acid
Flow Rate: 0.6 mL/min
Detection: RI @ 40 °C
Temperature: 60 °C
System: Shimadzu Prominence LC-20A System
Sample: 1. Dp4+ 5. Lactic Acid
 2. Dp3 6. Glycerol
 3. Maltose 7. Acetic Acid
 4. Glucose 8. Ethanol



Shorten GC Fuel Quality Testing

Zebtron™ ZB-Bioethanol GC column can shorten your quality test down to 5 minutes! (See p. 99).



Extend Column Lifetime

Protect the Rezex column from the intrusion of the metal ions by using Phenex™ Syringe Filters and SecurityGuard™. The filters and SecurityGuard guard cartridge system work by trapping metal ions, such as calcium, magnesium, and iron, which can damage the column and cause it to lose or change separation efficiency. (See pp. 12 and 251).

Rezex™ vs. Bio-Rad Aminex®

Phenomenex guarantees satisfaction when using Rezex HPLC columns. As illustrated below, Rezex offers advantages that enhance chromatographic results, increase throughput, and simplify quantitation.

Easier, Accurate Quantitation

Due to improved peak shape

Saccharides

Conditions for both columns:

Column: Rezex RCM-Monosaccharide
Aminex HPX-87C

Dimensions: 300 x 7.8 mm

Mobile Phase: Water


Flow Rate: 0.6 mL/min

Detection: ELSD

Temperature: 80 °C

Sample: 1. Melezitose 4. Mannose
2. Maltose 5. Fructose
3. Glucose 6. Ribitol

Superior
Peak Shape



Orange Juice

Conditions for both columns:

Column: Rezex RCM-Monosaccharide
Aminex HPX-87C

Dimensions: 300 x 7.8 mm


Mobile Phase: Water

Flow Rate: 0.6 mL/min

Detection: ELSD

Temperature: 80 °C

Sample: 1. Sucrose
2. Glucose
3. Fructose



Longer Column Lifetime

Due to lower backpressures

Aliphatic Acids

Conditions for both columns:

Column: Rezex ROA-Organic Acid
Aminex HPX-87H

Dimensions: 300 x 7.8 mm


Mobile Phase: 0.005 N H₂SO₄

Flow Rate: 0.5 mL/min

Detection: UV @ 210 nm

Temperature: 40 °C

Sample: 1. Oxalic Acid 4. Succinic Acid
2. Citric Acid 5. Formic Acid
3. Tartaric Acid 6. Acetic Acid



Baseline Separation of Critical Sample Components

Due to improved resolution

Sugars

Conditions for both columns:

Column: Rezex RCM-Monosaccharide
Aminex HPX-87C

Dimensions: 300 x 7.8 mm


Mobile Phase: Water

Flow Rate: 0.6 mL/min

Detection: ELSD

Temperature: 80 °C


Sample: 1. Sucrose 4. Fructose
2. Glucose 5. Sorbitol
3. Galactose



Comparative separations may not be representative of all applications.


Organic Acids

Column: Rezex RQA-Organic Acid
Dimensions: 300 x 7.8 mm
Part No.: 00H-0138-KO
Mobile Phase: 0.005 N Sulfuric Acid
Flow Rate: 0.5 mL/min
Detection: UV @ 210 nm
Temperature: 55 °C
Sample: 1. Oxalic
 2. Citric
 3. Tartaric
 4. Succinic
 5. Formic
 6. Acetic




Saccharides

Column: Rezex RCM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0130-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 85 °C
Sample: 1. Melezitose
 2. Maltose
 3. Glucose
 4. Mannose
 5. Fructose
 6. Ribitol




Apple Juice

Column: Rezex RCM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0130-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 75 °C
Sample: 1. Sucrose
 2. Glucose
 3. Fructose




Carbohydrates

Column: Rezex RKP-Potassium
Dimensions: 300 x 7.8 mm
Part No.: 00H-3252-KO
Mobile Phase: Water
Flow Rate: 0.4 mL/min
Detection: RI
Temperature: 85 °C
Sample: 1. Maltotriose
 2. Maltose
 3. Glucose




Oligosaccharides

Column: Rezex RSO-Oligosaccharide
Dimensions: 200 x 10 mm
Part No.: 00P-0133-NO
Mobile Phase: Water
Flow Rate: 0.3 mL/min
Detection: RI
Temperature: 75 °C
Sample: Malto-Oligosaccharides as shown




Saccharides

Column: Rezex RPM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0135-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 75 °C
Sample: 1. Stachyose
 2. Maltose
 3. Glucose
 4. Xylose
 5. Galactose
 6. Fructose
 7. Meso-Erythritol
 8. Mannitol
 9. Salicin
 10. Xylitol
 11. Sorbitol




Fermentation Broth

Column: Rezex RCM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0130-KO
Mobile Phase: Water
Flow Rate: 0.5 mL/min
Detection: ELSD
Temperature: 80 °C
Sample: 1. Sucrose
 2. Glucose
 3. Galactose
 4. Fructose




Mannitol and Sorbitol: USP Method (L19 Column)

Column: Rezex RPM-Monosaccharide
Dimensions: 100 x 7.8 mm
Part No.: 00D-0135-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: ELSD
Temperature: 80 °C
Sample: 1. Mannitol
 2. Sorbitol




Carbohydrates

Column: Rezex RNM-Carbohydrate
Dimensions: 300 x 7.8 mm
Part No.: 00H-0136-KO
Mobile Phase: Water
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: 75 °C
Sample: 1. Stachyose
 2. Cellobiose
 3. Glucose
 4. Fructose
 5. Arabinose
 6. Ribose
 7. Salicin
 (1% solution of each)




Carboxylic Acids

Column: Rezex ROA-Organic Acid
Dimensions: 300 x 7.8 mm
Part No.: 00H-0138-KO
Mobile Phase: Water + 0.5 % Trifluoroacetic Acid
Flow Rate: 1.0 mL/min
Detection: RI
Temperature: 40 °C
Sample: 1. Acetylene Carboxylic Acid
 2. Maleic Acid
 3. Succinic Acid
 4. Fumaric Acid




Food Softeners

Column: Rezex RCM-Monosaccharide
Dimensions: 300 x 7.8 mm
Part No.: 00H-0130-KO
Mobile Phase: Water
Flow Rate: 0.5 mL/min
Detection: RI
Temperature: 60 °C
Sample: 1. Glycerol
 2. Methoxypolyethylene Glycol
 3. Triethylene Glycol
 4. Sorbitol
 5. Urea



Amino Sugars

Column: Rezex ROA-Organic Acid
Dimensions: 300 x 7.8 mm
Part No.: 00H-0138-KO
Mobile Phase: 1 % Phosphoric Acid
Flow Rate: 0.6 mL/min
Detection: RI
Temperature: Ambient
Sample: 1. Glucose
 2. N-Acetylglucosamine
 3. N-Acetylgalactosamine



Specifications and Operating Recommendations

	RCM-Monosaccharide	RSO-Oligosaccharide	RNO-Oligosaccharide	RNM-Carbohydrate	RAM-Carbohydrate
Part Number	00H-0130-K0	00P-0133-N0	00P-0137-N0	00H-0136-K0	00H-0131-K0
Ionic Form	Calcium	Silver	Sodium	Sodium	Silver
Standard Dimensions	300 x 7.8 mm	200 x 10 mm	200 x 10 mm	300 x 7.8 mm	300 x 7.8 mm
Matrix	Sulfonated Styrene-divinylbenzene				
Cross Linking	8 %	4 %	4 %	8 %	8 %
Particle Size	8 µm	12 µm	12 µm	8 µm	8 µm
Min. Efficiency (p/m) based on last peak	35,000	N/A	N/A	30,000	35,000
Typical Pressure (psi @ Testing Flow Rate)	260	115	130	170	285
Max. Pressure (psi @ Max Flow Rate)	1,000	300	300	1,000	1,000
Max. Flow Rate (mL/min)	1.0 (see pressure)	0.3	0.3	1.0	1.0
Max. Temperature (°C)	85	85	85	85	85
Typical Mobile Phase	Water	Water	Water	Water	Water
pH Range	Neutral	Neutral	Neutral	Neutral	Neutral
Guard Column Part No.	03B-0130-K0	03R-0133-N0	03R-0137-N0	03B-0136-K0	03B-0131-K0

Cleaning, Regeneration and Storage

Organic Modifiers (Max)			5 % Methanol, IPA, EtOH		
Inorganic Modifiers	5 % CaSO ₄ , Ca(NO ₃) ₂ , CaCl ₂	5 % Silver Nitrate	5 % Sodium Salts	5 % Sodium Salts	2 % Silver Nitrate
Avoid 	Acids, Bases, Non-Calcium Salts/ Metal Ions, >30 % Acetonitrile	Acids, Bases, Non-Silver Salts/ Metal Ions, >30 % Acetonitrile	Acids, Bases, Non-Sodium Salts/ Metal Ions, >30 % Acetonitrile	Acids, Bases, Non-Sodium Salts/ Metal Ions, >30 % Acetonitrile	Acids, Bases, Non-Silver Salts/ Metal Ions, >30 % Acetonitrile
Cleaning Solvent	100 % Water	100 % Water	100 % Water	100 % Water	100 % Water
Flow Rate(mL/min)	0.4	0.1	0.1	0.4	0.4
Temperature (°C)	85	85	85	85	85
Duration (hrs)	12	12	12	12	12
Regeneration Solvent	0.1 M Ca(NO ₃) ₂	0.1 M AgNO ₃	0.1 M NaNO ₃	0.1 M NaNO ₃	0.1 M AgNO ₃
Flow Rate (mL/min)	0.2	0.1	0.2	0.2	0.2
Temperature (°C)	85	85	85	85	85
Duration (hrs)	4-16	4-16	4-16	4-16	4-16
Ship/Storage Solvent	Water	Water	Water	Water	Water

	RPM-Monosaccharide	RHM-Monosaccharide	ROA-Organic Acid	RFQ-Fast Acid	RCU-Sugar Alcohols
Part Number	00H-0135-K0	00H-0132-K0	00H-0138-K0	00D-0223-K0	00G-0130-D0
Ionic Form	Lead	Hydrogen	Hydrogen	Hydrogen	Calcium
Standard Dimensions	300 x 7.8 mm	300 x 7.8 mm	300 x 7.8 mm	100 x 7.8 mm	250 x 4.0 mm
Matrix	Sulfonated Styrene-divinylbenzene				
Cross Linking	8 %	8 %	8 %	8 %	8 %
Particle Size	8 µm	8 µm	8 µm	8 µm	8 µm
Min. Efficiency (p/m) based on last peak	35,000	35,000	50,000 (Acetic Acid)	30,000	12,000
Typical Pressure (psi @ Testing Flow Rate)	190	275	580	365	90
Max. Pressure (psi @ Max Flow Rate)	1,000	1,000	1,000	1,000	1,000
Max. Flow Rate (mL/min)	1.0	1.0	1.0	1.0	0.5
Max. Temperature (°C)	85	85	85	85	85
Typical Mobile Phase	Water	Water	0.005 N H ₂ SO ₄	0.005 N H ₂ SO ₄	Water
pH Range	Neutral	1-8	1-8	1-8	Neutral
Guard Column Part No.	03B-0135-K0	03B-0132-K0	03B-0138-K0	03B-0223-K0	03A-0130-D0

Cleaning, Regeneration and Storage

Organic Modifiers (Max)			5 % Methanol, IPA, EtOH		
Inorganic Modifiers	5 % Lead Nitrate	5 % HNO ₃ , H ₃ PO ₄	5 % HNO ₃ , H ₃ PO ₄	5 % HNO ₃ , H ₃ PO ₄	5 % CaSO ₄ , Ca(NO ₃) ₂ , CaCl ₂
Avoid 	Acids, Bases, Non-Lead Salts/ Metal Ions, >30 % Acetonitrile	Acids, Bases, Salts/ Metal Ions, >30 % Acetonitrile	Acids, Bases, Salts, Metal Ions, pH > 3, >30 % Acetonitrile	Acids, Bases, Salts, Metal Ions, pH > 3, >30 % Acetonitrile	Acids, Bases, Non-Calcium Salts, or Metal Ions, >30 % Acetonitrile
Cleaning Solvent	100 % Water	100 % Water	100 % Water	100 % Water	100 % Water
Flow Rate(mL/min)	0.4	0.4	0.4	0.4	0.1
Temperature (°C)	85	85	85	85	85
Duration (hrs)	12	12	12	12	12
Regeneration Solvent	0.1 M Pb(NO ₃) ₂	0.025 M H ₂ SO ₄	0.025 M H ₂ SO ₄	0.025 M H ₂ SO ₄	0.1 M Ca(NO ₃) ₂
Flow Rate (mL/min)	0.2	0.2	0.2	0.2	0.1
Temperature (°C)	85	85	85	85	85
Duration (hrs)	4-16	4-16	4-16	4-16	4-16
Ship/Storage Solvent	Water	Water	0.005 N H ₂ SO ₄	0.005 N H ₂ SO ₄	Water

Retention Times for Some Carbohydrates and Sugar Alcohols

Counter Ion	Analyte	RAM Ag ⁺	RCM Ca ⁺²	RNM Na ⁺	RHM H ⁺	RPM Pb ⁺²
Adonitol (Ribitol)		11.54	14.93	11.10	11.11	20.15
D-Altrose		11.95	12.71	11.45	10.21	15.82
D-(-)-Arabinose		13.01	13.56	12.65	11.24	16.47
D-(+)-Cellobiose		8.86	8.60	8.49	8.02	11.00
D-(+)-Digitoxose		11.90	13.82	11.39	12.59	15.32
Dulcitol		11.64	21.61	11.10	10.71	33.25
Meso-Erythritol		12.31	15.49	11.78	12.14	19.82
D-(-) Fructose		12.05	13.65	11.76	10.31	17.71
L-(-)-Fucose		12.75	13.19	12.30	11.65	16.19
D-(+)-Galactose		11.87	11.73	11.47	10.19	14.94
Gentiobiose		8.70	8.40	8.40	7.87	10.53
D-(+)-Glucose		11.04	10.37	10.71	9.62	12.92
Inositol		12.59	13.35	12.14	9.98	18.87
Isomaltose		9.11	8.74	8.76	8.02	11.28
Lactose		9.27	9.03	8.78	8.32	11.89
Lactulose		9.75	10.32	9.23	8.57	13.95
D- Lyxose		12.41	14.06	11.98	10.68	16.66
D- Maltose		9.16	8.81	8.75	8.18	11.59
Maltotriose		8.27	8.10	7.94	7.51	11.02
Maltulose		9.25	9.47	8.82	8.27	12.40
D- Mannitol		11.36	17.82	10.80	10.59	24.90
D-(+)-Mannose		12.04	12.04	11.54	10.16	16.39
Melibiose		9.26	9.04	8.82	8.14	11.97
D-(+)-Melezitose		8.00	7.93	7.66	7.54*	9.94
D-(+)-Raffinose		8.10	8.16	7.76	7.88*	10.28
L-(+)-Rhamnose		11.50	12.18	11.00	10.90	14.47
D-(-)-Ribose		14.59	23.38	14.34	11.42	33.48
Salicin		18.51	18.58	17.36	14.98	26.81
D-Sorbitol		11.91	22.45	11.39	10.83	35.97
Stachyose		7.60	7.59	7.30	7.27	9.72
Sucrose		9.03	8.71	8.65	9.24*	11.00
Trehalose		8.91	8.72	8.49	8.32	11.01
Xylitol		12.69	22.01	12.16	11.78	32.38
D-(+)-Xylose		12.06	11.62	11.68	10.24	13.84

* Partial hydrolysis results.

Conditions:

Dimensions: 300 x 7.8 mm
Mobile Phase: Water (degassed)
Flow Rate: 0.6 mL/min
Temperature: 80 °C
Detection: RI @ 40 °C

If you are not completely satisfied with the performance of any Rezex column, as compared to a competing product of the same size and phase, simply return the Rezex column with your comparative data within 45 days for a FULL REFUND.

Column Cross Reference Chart

Phenomenex Rezex™	Bio-Rad® Aminex®	Supelco® SUPELCOGEL™	Waters® Sugar-Pak™
RCM-Monosaccharide	HPX-87C 125-0095	Supelcogel Ca	Sugar-Pak 1
RHM-Monosaccharide	HPX-87H 125-0140	Supelcogel C-610H & H	N/A
RPM-Monosaccharide	HPX-87P 125-0098	Supelcogel Pb	N/A
RNM-Carbohydrate	HPX-87N 125-0143	N/A	N/A
RSO-Oligosaccharide	HPX-42A 125-0097	Supelcogel Ag1 & Ag2	N/A
ROA-Organic Acid	HPX-87H 125-0140	Supelcogel C-610H & H	N/A
RFQ-Fast Acid	Fast Acid 125-0100	N/A	N/A
RKP-Potassium	HPX-87K 125-0142	Supelcogel K	N/A
RCU-USP Sugar Alcohols	Sugar Alcohols 125-0094	N/A	N/A

Ordering Information

Columns						Guards			SecurityGuard™ Cartridges (mm)
Description	Part No.	Cross Linkage	Ionic Form	Size (mm)	Price	Part No.	Size (mm)	Price	4 x 3.0* /10pk
RCM-Monosaccharide	00F-0130-K0	8 %	Calcium	150 x 7.8		03B-0130-K0	50 x 7.8		AJ0-4493
RCM-Monosaccharide	00H-0130-K0	8 %	Calcium	300 x 7.8		03B-0130-K0	50 x 7.8		AJ0-4493
RHM-Monosaccharide	00H-0132-K0	8 %	Hydrogen	300 x 7.8		03B-0132-K0	50 x 7.8		AJ0-4490
RAM-Carbohydrate	00H-0131-K0	8 %	Silver	300 x 7.8		—	—	—	AJ0-4491
RSO-Oligosaccharide	00P-0133-N0	4 %	Silver	200 x 10.0		03R-0133-N0	60 x 10.0		—
RNO-Oligosaccharide	00P-0137-N0	4 %	Sodium	200 x 10.0		03R-0137-N0	60 x 10.0		—
RPM-Monosaccharide (for USP procedure)	00H-0135-K0	8 %	Lead	300 x 7.8		03B-0135-K0	50 x 7.8		AJ0-4492
	00D-0135-K0	8 %	Lead	100 x 7.8		03B-0135-K0	50 x 7.8		AJ0-4492
RNM-Carbohydrate	00H-0136-K0	8 %	Sodium	300 x 7.8		03B-0136-K0	50 x 7.8		—
ROA-Organic Acid	00F-0138-E0	8 %	Hydrogen	150 x 4.6		—	—	—	AJ0-4490
ROA-Organic Acid	00G-0138-E0	8 %	Hydrogen	250 x 4.6		—	—	—	AJ0-4490
ROA-Organic Acid	00F-0138-K0	8 %	Hydrogen	150 x 7.8		03B-0138-K0	50 x 7.8		AJ0-4490
ROA-Organic Acid	00H-0138-K0	8 %	Hydrogen	300 x 7.8		03B-0138-K0	50 x 7.8		AJ0-4490
RKP-Potassium	00H-3252-K0	8 %	Potassium	300 x 7.8		—	—	—	—
RFQ-Fast Acid	00D-0223-K0	8 %	Hydrogen	100 x 7.8		03B-0223-K0	50 x 7.8		AJ0-4490
RCU-USP Sugar Alcohols	00G-0130-D0	8 %	Calcium	250 x 4.0		03A-0130-D0	30 x 4.0		AJ0-4493

for ID: 3.2-8.0 mm

*SecurityGuard Analytical Cartridges require universal holder Part No.: KJ0-4282



For Column Heaters, see p. 366



For our full line of Column Performance Check Standards, see pp. 372-373