



# Zebtron™ Inferno™

Rugged High Temperature GC Columns

ZB-1HT

ZB-5HT

**NEW!** ZB-35HT  
**NEW!** ZB-XLB-HT



 **phenomenex**<sup>®</sup>  
...breaking with tradition<sup>SM</sup>



# Manufacturing Award Winning Solutions for GC Analysis




In the field of gas chromatography, not much has changed in the last 15 years. However, technology is our passion and Phenomenex is gaining international recognitions for its GC products that are specially designed for today's most challenging methods.

For chemists struggling with high boilers, contaminants, or carry-overs, Phenomenex developed the Zebtron™ Inferno™ line. This line of GC columns have revolutionized GC analysis by allowing chemists to run their analysis at oven temperatures that were never before attainable - all without having to use a metal column!

The Zebtron ZB-1HT and ZB-5HT Inferno columns are rugged, mid-polarity columns that allow chemists to push their standard analysis up to 430 °C for increased lab productivity and improved performance. The newly developed Zebtron ZB-35HT and ZB-XLB-HT Inferno columns provide high temperature capability up to 400 °C for those using more polar columns.

guarantee


If Zebtron does not provide you with equivalent separations as compared to any other GC column of the same phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.





# Zebtron<sup>TM</sup> Inferno<sup>TM</sup>

Rugged High Temperature GC Analysis



|    |                                                              |
|----|--------------------------------------------------------------|
| 4  | Challenges with Traditional GC Columns                       |
| 5  | Advanced Solutions for High Temperature Analysis             |
| 6  | Low Bleed                                                    |
| 7  | Increased Column Lifetime                                    |
| 8  | Ideal for High Boiling Compounds                             |
| 10 | Bake-off Contaminants and Eliminate Carry-overs!             |
| 11 | See What You Have Been Missing!                              |
| 12 | Widely Used Phases (Zebtron ZB-1HT and ZB-5HT)               |
| 13 | Mid-Polar Phases (Zebtron ZB-35HT and ZB-XLB-HT) <b>NEW!</b> |
| 15 | Applications                                                 |
|    | Environmental ..... 16                                       |
|    | Petrochemical ..... 21                                       |
|    | Biodiesel ..... 22                                           |
|    | Food and Beverage ... 23                                     |
|    | Chemical ..... 24                                            |
|    | Pharmaceutical ..... 25                                      |
| 25 | Ordering Information                                         |
|    | ZB-1HT ..... 12                                              |
|    | ZB-5HT ..... 12                                              |
|    | ZB-35HT ..... 14 <b>NEW!</b>                                 |
|    | ZB-XLB-HT ..... 14 <b>NEW!</b>                               |
|    | GC Accessories ..... 25                                      |

# Challenges with Traditional GC Columns

## Standard fused silica columns - Shorter lifetime

Unlike the Zebron™ Inferno™ GC columns, most fused silica columns are not developed to handle prolonged exposure to temperatures above 360 °C.

**At high temperatures, fused silica GC columns:**

- Become inflexible and brittle
- Break easily and are costly to replace
- Require more system downtime to replace columns
- Can't separate high molecular weight compounds

Deterioration of a standard column at high temperatures




## Metal columns – Difficult to use

Though metal columns have higher temperature limits than capillary columns, they are more difficult to use.


**Metal columns:**

- Inflexible and hard to use
- Require special tubing cutters for installation
- Prone to leaks from expansion and contraction
- Have active surface, giving poor peak shapes for acids and bases
- Incompatible with MS detectors



# Advanced Solutions for High Temperature Analysis

Zebtron™ Inferno™ columns offer enhanced performance where other columns fail. The two critical factors that allow for this are:



## 1. Advanced ESC™ Bonding Technology

At high temperature ranges, the stability of standard GC columns will deteriorate, resulting in increased bleed. Zebtron's Engineered Self Cross-linking™ (ESC) bonding technology reinforces the stationary phase for enhanced column durability and extremely low bleed levels. The result - our Inferno columns allow for the flexibility to perform GC analysis at high temperatures while providing low bleed and longer column lifetime.


## 2. Special High Temperature Polyimide Coating

Standard polyimide resin pyrolyzes at temperatures above 360 °C, making the tubing unstable. The Zebtron Inferno columns utilize a special dark polyimide resin that shows minimal thermal degradation even at temperatures up to 430 °C.\* This results in longer column lifetime at elevated temperatures.

| Phase     | Temperature Limit |        |
|-----------|-------------------|--------|
| ZB-1HT    | 430 °C            |        |
| ZB-5HT    | 430 °C            |        |
| ZB-35HT   | <b>NEW!</b>       | 400 °C |
| ZB-XLB-HT | <b>NEW!</b>       | 400 °C |

## Flexible Even at 430 °C!\*\*

The Zebtron Inferno columns remain flexible and easy to work with, even after being exposed to many hours at extreme temperatures. Its high temperature stability offers the flexibility of a non-metal column while providing long lifetime, low bleed, and low activity.



\* Zebtron ZB-1HT and ZB-5HT Inferno columns have an upper temperature limit of 430 °C. Zebtron ZB-35HT and ZB-XLB-HT Inferno columns have an upper temperature limit of 400 °C.

\*\* Evaluated by performing 185 programmed temperature runs, total 23 hours at 430 °C. Polyimide tubing was still flexible as shown here.

# Extremely Low Bleed for Rugged Performance

- Longer column lifetime
- Better reproducibility
- Maximum sensitivity
- Accurate quantitation

## What is column bleed?

Column bleed is the loss of stationary phase as a result of impurities in the starting polymer or the decomposition of the phase at elevated temperatures.


## Why is column bleed undesirable?

Column bleed can be an indicator of the stability and lifetime a GC column will offer. When using sensitive detectors such as MS, bleed can impact method detection limits (MDLs) causing difficulties at the low calibration ranges.

## What sets Zebron's bleed apart from others?

Our proprietary ESC™ bonding technology and special polyimide tubing provides enhanced column durability for extremely low bleed levels.

## Zebron™ Inferno™ Columns Have Lower Bleed than Other Manufacturers†



† All columns were new, never used, prior to this testing. Columns were purchased either directly from the original manufacturer or through an authorized distributor. All testing was carefully controlled to ensure conditions were similar for all columns involved. The comparative data is not representative of every application.


# Zebtron™ Inferno™'s High Durability Provides Increased Column Lifetime

- Less instrument downtime
- More productivity
- Lower laboratory costs

## Standard GC Columns Cannot Withstand High Temperature

After prolonged exposure to temperatures above 360 °C, most traditional GC columns become brittle and inflexible, often spontaneously breaking in the middle of a run. Work is stopped in order to change out the column and perform system maintenance. Unlike these fragile columns, the Zebtron Inferno columns have a special polyimide coating that provides longer column lifetime at temperatures up to 430 °C. As a result, you will experience longer column lifetime, less instrument downtime, and higher productivity.

## Zebtron Inferno Columns Win In The Lifetime Test



\*All columns were new, never used, prior to this testing. Columns were purchased either directly from the original manufacturer or through an authorized distributor. All testing was carefully controlled to ensure conditions were similar for all columns involved. The comparative data is not representative of every application.

## How does the lifetime test work?

Hydrocarbons are a good way to measure the stability and lifetime of a non-polar column. Because its interaction with the phase is mostly based on London Dispersion forces, any change in retention time is correlated with phase loss. This will result in increased bleed and poor reproducibility.

In the lifetime study above, the Zebtron ZB-5HT Inferno column has twice the lifetime as the two comparative columns. For the test, all columns were held at 400 °C for 2 hours and then the oven was lowered to 120 °C for Pentadecane analysis. The Varian® HT column broke just after 40 hours at 400 °C. The Agilent® J&W® DB-5ht column had the same retention for Pentadecane at 40 hours as the ZB-5HT at 100 hours.




# Greater Performance for High Boiling Compounds

- Better peak shape
- Improved separation
- More analytes detected

Zebron™ Inferno™ GC columns have the inertness and temperature stability needed for separating high molecular weight compounds. These high boiling compounds require a thermally rugged column with low bleed levels. Zebron Inferno columns have the high temperature stability and durability to provide optimal separation of these compounds.

## Great Separation of High Boiling Hydrocarbons (ASTM Method D6352)



### Easily Elute Hydrocarbons Up To C90

Zebron's ZB-1HT Inferno can operate at 430 °C, allowing C90 to elute on temperature ramp to meet ASTM requirements


**Column:** Zebron ZB-1HT Inferno  
**Dimensions:** 5 meter x 0.53 mm x 0.10 µm  
**Part No.:** 7AK-G014-02  
**Injection:** On-Column @ 43 °C, 0.1 µL  
**Carrier Gas:** Helium @ 4.4 mL/min (constant flow)  
**Oven Program:** 40 °C for 0.5 min to 430 °C @ 20 °C/min for 10 min  
**Detector:** FID @ 430 °C

Note: Sample was a combination of PolyWax 655 and retention time markers C8-C40 in CS2/Chloroform

| Sample:     |             |         |         |         |         |
|-------------|-------------|---------|---------|---------|---------|
| 1. C10      | 11. Phytane | 21. C28 | 31. C38 | 41. C56 | 51. C76 |
| 2. C11      | 12. C19     | 22. C29 | 32. C39 | 42. C58 | 52. C78 |
| 3. C12      | 13. C20     | 23. C30 | 33. C40 | 43. C60 | 53. C80 |
| 4. C13      | 14. C21     | 24. C31 | 34. C42 | 44. C62 | 54. C82 |
| 5. C14      | 15. C22     | 25. C32 | 35. C44 | 45. C64 | 55. C84 |
| 6. C15      | 16. C23     | 26. C33 | 36. C46 | 46. C66 | 56. C86 |
| 7. C16      | 17. C24     | 27. C34 | 37. C48 | 47. C68 | 57. C88 |
| 8. C17      | 18. C25     | 28. C35 | 38. C50 | 48. C70 | 58. C90 |
| 9. Pristane | 19. C26     | 29. C36 | 39. C52 | 49. C72 |         |
| 10. C18     | 20. C27     | 30. C37 | 40. C54 | 50. C74 |         |



## Separation of Sterols from Margarine on a Zebron™ ZB-5HT Inferno™



App ID 16530

**Column:** Zebron ZB-5HT Inferno  
**Dimensions:** 30 meter x 0.25 mm x 0.10 µm  
**Part No.:** 7HG-G015-02  
**Injection:** Splitless @ 275 °C, 0.5 µL  
**Carrier Gas:** Helium @ 1.5 mL/min (constant flow)  
**Oven Program:** 220 °C to 350 °C @ 20 °C/min  
**Detector:** MSD @ 275 °C

**Sample:** 1. Campesterol  
2. Stigmasterol  
3. β-Sitosterol  
4. Betulin (IS)




# Bake-off Contaminants and Eliminate Carry-overs!

- Achieve better reproducibility
- Extend column lifetime
- Attain sharper peak shapes

## Eliminate Problems Associated with Dirty Samples

Dirty samples can linger on a GC column, eventually changing the selectivity and shortening the lifetime of a column. Since Zebron™ Inferno™ columns have extremely high temperature limits, it allows for the use of more aggressive column baking steps to remove persistent contaminants, thus prolonging column lifetimes.

Hydrocarbon contaminants are not detected by ECD!




Conditions for App ID 16306:

**Column:** Zebron ZB-5HT Inferno  
**Dimensions:** 30 meter x 0.32 mm x 0.10 µm  
**Part No.:** 7HM-G015-02  
**Injection:** Split 52:1 @ 330 °C, 1 µL  
**Carrier Gas:** Helium @ 1.9 mL/min (constant flow)  
**Oven Program:** 150 °C for 2 min to 420 °C @ 20 °C/min for 4.5 min  
**Detector:** ECD @ 380 °C  
**Sample:** 1. TCMX  
2. DCOB  
3. Aroclor 1016  
4. Mineral Oil

Note: Aroclor 1016 @ 100 ppb; Mineral Oil @ 1 % in isooctane

Hydrocarbon contaminants revealed using FID starting at 5 minutes



Conditions for App ID 18778:

**Column:** Zebron ZB-5HT  
**Dimensions:** 30 meter x 0.32 mm x 0.10 µm  
**Part No.:** 7HM-G015-02  
**Injection:** Split 5:1 @ 330 °C, 1 µL  
**Carrier Gas:** Helium @ 1.9 mL/min (constant flow)  
**Oven Program:** 150 °C for 2 min to 420 °C @ 20 °C/min for 4.5 min  
**Detector:** FID @ 435 °C  
**Sample:** Mineral Oil


Note: Mineral Oil is 1 % in isooctane

# See What You have been Missing!


- Uncover compounds that were previously undetected
- Elute compounds at higher temperatures
- Achieve more accurate, reproducible results

Many high boiling compounds will stay on a column and go undetected run after run. As a consequence, results are inaccurate because of missing data. Running at higher temperatures will allow these formerly undetected compounds to elute and be detected, giving more accurate and reproducible results.

## New Zebron™ ZB-35HT Inferno™ Column



## Restek® Rtx®-35 GC Column



# Available in Widely Used Phases

## ZB-1HT Inferno™

- Temperature Limits: -60 to 400/430 °C (Isothermal/TPGC)\*
- First non-metal columns stable to 430 °C
- Provides true boiling point separation for hydrocarbon distillation methods
- Individually tested for low bleed, MS certified

Alternative to any 100 % Dimethylpolysiloxane High-Temperature Phase:

• DB-1ht • MXT-1 SimDist • Petrocol 2887 • CP-SimDist

### Column Profile



### Applications

- High Boiling Petroleum Products
- Simulated Distillation Methods
- Long-chained Hydrocarbons
- High Molecular Weight Waxes
- Polymers/Plastics
- Motor Oils
- Diesel Fuel


guarantee

If Zebtron does not provide you with equivalent separations as compared to any other GC column of the same phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

### Ordering Information

#### Zebtron ZB-1HT Inferno GC Columns

| ID(mm)          | df(μm) | Temp. Limits °C | Part No.    | Price |
|-----------------|--------|-----------------|-------------|-------|
| <b>5-Meter</b>  |        |                 |             |       |
| 0.53            | 0.10   | -60 to 400/430  | 7AK-G014-02 |       |
| <b>15-Meter</b> |        |                 |             |       |
| 0.25            | 0.10   | -60 to 400/430  | 7EG-G014-02 |       |
| 0.32            | 0.10   | -60 to 400/430  | 7EM-G014-02 |       |
| 0.32            | 0.25   | -60 to 400/430  | 7EM-G014-11 |       |
| 0.53            | 0.15   | -60 to 400      | 7EK-G014-05 |       |
| <b>20-Meter</b> |        |                 |             |       |
| 0.18            | 0.18   | -60 to 400/430  | 7FD-G014-08 |       |
| <b>30-Meter</b> |        |                 |             |       |
| 0.25            | 0.10   | -60 to 400/430  | 7HG-G014-02 |       |
| 0.25            | 0.25   | -60 to 400/430  | 7HG-G014-11 |       |
| 0.32            | 0.10   | -60 to 400/430  | 7HM-G014-02 |       |
| 0.32            | 0.25   | -60 to 400/430  | 7HM-G014-11 |       |
| 0.53            | 0.15   | -60 to 400      | 7HK-G014-05 |       |



**Extend column lifetime. Add a Z-guard to your next Zebtron GC order.**

If you need a 5 in. cage, simply add a (-B) after the part number, e.g., 7HG-G014-11-B. Some exceptions may apply.  
\*0.53 mm ID columns are rated to 400 °C max operational temperature.

## ZB-5HT Inferno™

- Temperature Limits: -60 to 400/430 °C (Isothermal/TPGC)\*
- First non-metal columns stable to 430 °C
- Individually tested for low bleed, MS certified

Alternative to Any 5 %-Phenyl-95 % Dimethylpolysiloxane High Temperature Phase:

• DB-5ht • VF-5ht • HT-5 • Stx-5HT • XTI-5HT

### Column Profile




### Applications

- High Boiling Petroleum Products
- Triglycerides
- Simulated Distillation Methods
- Diesel Fuel
- Long-chained Hydrocarbons
- Motor Oils
- Polymers/Plastics
- Surfactants
- High Molecular Weight Waxes

### Ordering Information

#### Zebtron ZB-5HT Inferno GC Columns

| ID(mm)                                                  | df(μm) | Temp. Limits °C | Part No.        | Price |
|---------------------------------------------------------|--------|-----------------|-----------------|-------|
| <b>15-Meter</b>                                         |        |                 |                 |       |
| 0.25                                                    | 0.10   | -60 to 400/430  | 7EG-G015-02     |       |
| 0.25                                                    | 0.25   | -60 to 400/430  | 7EG-G015-11     |       |
| 0.32                                                    | 0.10   | -60 to 400/430  | 7EM-G015-02     |       |
| 0.32                                                    | 0.25   | -60 to 400/430  | 7EM-G015-11     |       |
| 0.53                                                    | 0.15   | -60 to 400      | 7EK-G015-05     |       |
| <b>15-Meter with 2-Meter Spliced Guard (0.53 mm ID)</b> |        |                 |                 |       |
| 0.32                                                    | 0.10   | -60 to 400/430  | 7EM-G015-02-GST |       |
| <b>20-Meter</b>                                         |        |                 |                 |       |
| 0.18                                                    | 0.18   | -60 to 400/430  | 7FD-G015-08     |       |
| <b>30-Meter</b>                                         |        |                 |                 |       |
| 0.25                                                    | 0.10   | -60 to 400/430  | 7HG-G015-02     |       |
| 0.25                                                    | 0.25   | -60 to 400/430  | 7HG-G015-11     |       |
| 0.32                                                    | 0.10   | -60 to 400/430  | 7HM-G015-02     |       |
| 0.32                                                    | 0.25   | -60 to 400/430  | 7HM-G015-11     |       |
| 0.53                                                    | 0.15   | -60 to 400      | 7HK-G015-05     |       |



**Extend column lifetime. Add a Z-guard to your next Zebtron GC order.**


If you need a 5 in. cage, simply add a (-B) after the part number, e.g., 7HG-G015-11-B. Some exceptions may apply.  
\*0.53 mm ID columns are rated to 400 °C max operational temperature.



# Two New Innovative Mid-Polar Phases:


## Zebron™ ZB-35HT and ZB-XLB-HT Inferno™

Unlike standard mid-polar GC columns that fail at high oven temperatures, the NEW Zebron ZB-35HT and ZB-XLB-HT Inferno columns have the durability and ruggedness to perform well at extreme temperatures. With an astonishing upper temperature limit of 400 °C, both columns will provide the separation of high boiling compounds at extremely low bleed levels.



## Zebron Inferno Provides the Lowest Bleed†

**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Injection:** Split 20:1 @ 200 °C, 1 µL  
**Carrier Gas:** Helium @ 1.7 mL/min (constant flow)  
**Oven Program:** 100 °C to 320 °C @ 30 °C/min for 5 min to 340 °C @ 20 °C/min for 5 min to 360 °C @ 20 °C/min for 5 min to 380 °C @ 20 °C/min for 5 min to 400 °C @ 20 °C/min for 5 min to 100 °C @ 30 °C/min for 8 min  
**Detector:** FID @ 405 °C



†All columns were new, never used, prior to this testing. Columns were purchased either directly from the original manufacturer or through an authorized distributor. All testing was carefully controlled to ensure conditions were similar for all columns involved. The comparative data is not representative of every application.

**NEW!**

## Zebtron ZB-35HT Inferno

- First non-metal, 35 %-Phenyl columns stable to 400 °C
- Great for high molecular weight compounds
- Individually tested for low bleed, MS certified
- Equivalent to USP Phase G42
- Temperature Limits: -60 to 400 °C (Isothermal/TPGC)


**guarantee**

If Zebtron does not provide you with equivalent separations as compared to any other GC column of the same phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Alternative to Any 35 %-Phenyl-65 % Dimethylpolysiloxane High Temperature Phase:

- |           |         |            |           |          |            |         |          |
|-----------|---------|------------|-----------|----------|------------|---------|----------|
| • DB-35ms | • AT-35 | • SPB-608  | • HP-35ms | • MDN-35 | • Sup-Herb | • DB-35 | • Rtx-35 |
| • SPB-35  | • OV-11 | • Rtx-35MS | • HP-35   | • 007-11 | • BPX608   | • EC-35 | • BPX35  |

### Column Profile



### Ordering Information

#### Zebtron ZB-35HT Inferno GC Columns

| ID(mm)          | df(μm) | Temp. Limits °C | Part No.    | Price |
|-----------------|--------|-----------------|-------------|-------|
| <b>15-Meter</b> |        |                 |             |       |
| 0.25            | 0.25   | -60 to 400      | 7EG-G025-11 |       |
| <b>30-Meter</b> |        |                 |             |       |
| 0.25            | 0.25   | -60 to 400      | 7HG-G025-11 |       |
| 0.32            | 0.25   | -60 to 400      | 7HM-G025-11 |       |

### Applications

- |                                          |                   |
|------------------------------------------|-------------------|
| • Aroclors                               | • Pesticides      |
| • Semi-volatiles                         | • Drugs of Abuse  |
| • Amines                                 | • Pharmaceuticals |
| • EPA Methods 508, 608, 8081, 8141, 8151 | • Steroids        |
|                                          | • Chemicals       |



**ZB-35HT Test Mix** Part No.: AG0-5156

**NEW!**

## Zebtron ZB-XLB-HT Inferno

- Non-metal si-arylene GC column stable to 400 °C
- Good tool for general screening to identify unknown samples
- Individually tested for low bleed, MS certified
- Temperature Limits: -60 to 400 °C (Isothermal/TPGC)



Similar\* to:

- |          |          |           |
|----------|----------|-----------|
| • MDN-12 | • DB-XLB | • Rtx-XLB |
|----------|----------|-----------|

\*not exact equivalent, selectivity might be different

### Column Profile



### Ordering Information

#### Zebtron ZB-XLB-HT Inferno GC Columns

| ID(mm)          | df(μm) | Temp. Limits °C | Part No.    | Price |
|-----------------|--------|-----------------|-------------|-------|
| <b>15-Meter</b> |        |                 |             |       |
| 0.25            | 0.10   | -60 to 400      | 7EG-G024-02 |       |
| 0.25            | 0.25   | -60 to 400      | 7EG-G024-11 |       |
| 0.32            | 0.10   | -60 to 400      | 7EM-G024-02 |       |
| <b>30-Meter</b> |        |                 |             |       |
| 0.25            | 0.10   | -60 to 400      | 7HG-G024-02 |       |
| 0.25            | 0.25   | -60 to 400      | 7HG-G024-11 |       |
| 0.32            | 0.25   | -60 to 400      | 7HM-G024-11 |       |

### Applications

- |                                    |               |
|------------------------------------|---------------|
| • Polychlorinated Biphenyls (PCBs) | • Herbicides  |
| • Pesticides                       | • EPA Methods |



**ZB-XLB-HT Test Mix** Part No.: AG0-7578

The background of the entire page is a composite image. The top half features laboratory glassware, including a beaker, a graduated cylinder, and a test tube, some containing liquids. The bottom half shows a close-up of green, leafy plants. The entire image is overlaid with a semi-transparent orange gradient.

# Suited for a Wide Range of Applications

## **For Many Industries:**

Environmental .....pp.16-20  
Petrochemical .....p.21  
Biodiesel .....p.22  
Food/Beverage .....p.23  
Chemicals .....p.24  
Pharmaceuticals.....p.24

## **Dual Column Analysis**

EPA 552.2..... p.16  
EPA 8081 ..... p.18  
EPA 8151 ..... p.18  
EPA 8270..... p.17

## **High Molecular Weight Compounds**

Heavy PAHs..... p.19  
Long chained hydrocarbons..... p.21  
Oils..... pp.21, 23  
Triglycerides..... p.23

## **Dirty Samples**



Oils..... pp.21, 23

## **Low Level Analysis**

Environmental contaminants .....pp.16-20  
Pesticides .....pp.16-20  
Biodiesel .....p.22  
Pharmaceuticals.....p.24

# Environmental Applications

## EPA Method 552.2 - Haloacetic Acids



Conditions for both columns:

**Injection:** Split 5:1 @ 250 °C, 1 µL

**Carrier Gas:** Helium @ 4.0 mL/min (constant flow)

**Oven Program:** 40 °C for 1 min to 55 °C @ 5 °C/min for 5 min to 70 °C @ 10 °C/min for 9 min to 250 °C @ 30 °C/min for 2 min


**Detector:** ECD @ 300 °C

Note: Analytes are methylated and 100 ppb in isooctane.

**Sample:** 1. Chloroacetic Acid  
2. Bromoacetic Acid  
3. Dichloroacetic Acid  
4. Dalapon  
5. Trichloroacetic Acid  
6. 1,2,3-Trichloropropane

7. Bromochloroacetic Acid  
8. Bromodichloroacetic Acid  
9. Dibromoacetic Acid  
10. Chlorodibromoacetic Acid  
11. 2,3-Dibromopropionate  
12. Tribromoacetic Acid

## DIN Method 51527: Aroclors 1242



**Column:** Zebron ZB-XLB-HT Inferno™

**Dimensions:** 30 meter x 0.25 mm x 0.25 µm

**Part No.:** 7HG-G024-11

**Injection:** Split 2:1 @ 250 °C, 1 µL

**Carrier Gas:** Helium @ 1.5 mL/min (constant flow)

**Oven Program:** 50 °C for 0.5 min to 210 °C @ 40 °C/min for 3 min to 230 °C @ 30 °C/min for 5 min to 250 °C @ 30 °C/min for 5 min to 320 °C @ 40 °C/min for 2 min


**Detector:** ECD @ 350 °C

Note: Total concentration of aroclors was 90 ppm in isooctane.

**Sample:**

|            |             |             |
|------------|-------------|-------------|
| 1. TCMX    | 8. BZ# 118  | 15. BZ# 156 |
| 2. BZ# 31  | 9. BZ# 153  | 16. BZ# 180 |
| 3. BZ# 28  | 10. BZ# 114 | 17. BZ# 157 |
| 4. BZ# 52  | 11. BZ# 105 | 18. BZ# 169 |
| 5. BZ#101  | 12. BZ# 138 | 19. BZ# 169 |
| 6. BZ# 77  | 13. BZ# 126 | 20. BZ# 189 |
| 7. BZ# 123 | 14. BZ# 167 | 21. DCB     |

## DIN Method 51527: Aroclors 1242



**Column:** Zebron ZB-35HT Inferno™

**Dimensions:** 30 meter x 0.25 mm x 0.25 µm

**Part No.:** 7HG-G025-11

**Injection:** Splitless @ 250 °C, 1 µL

**Carrier Gas:** Helium @ 0.5 mL/min (constant flow)


**Oven Program:** 50 °C for 0.5 min to 220 °C @ 30 °C/min to 320 °C @ 6 °C/min for 2 min

**Detector:** ECD @ 350 °C

**Sample:**

|            |             |             |
|------------|-------------|-------------|
| 1. TCMX    | 8. BZ# 118  | 15. BZ# 156 |
| 2. BZ# 31  | 9. BZ# 153  | 16. BZ# 180 |
| 3. BZ# 28  | 10. BZ# 114 | 17. BZ# 157 |
| 4. BZ# 52  | 11. BZ# 105 | 18. BZ# 169 |
| 5. BZ# 101 | 12. BZ# 138 | 19. BZ# 170 |
| 6. BZ# 77  | 13. BZ# 126 | 20. BZ# 189 |
| 7. BZ# 123 | 14. BZ# 167 | 21. DCB     |

## EPA Method 8270 - Semi-volatile Organics



App ID 18785


**Column:** Zebron ZB-XLB-HT Inferno™  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G024-11  
**Injection:** Pulsed @ 280 °C, 1 µL  
**Carrier Gas:** Helium @ 1 mL/min (constant flow)  
**Oven Program:** 40 °C for 1.5 min to 260 °C @ 20 °C/min to 290 °C @ 5 °C/min to 350 °C @ 25 °C/min for 2 min  
**Detector:** MSD @ 350 °C; 45-450 amu  
 Note: Analytes are 20 ppm in dichloromethane

|                                 |                                |                                 |                                     |
|---------------------------------|--------------------------------|---------------------------------|-------------------------------------|
| <b>Sample:</b> 1. Pyridine      | 36. 2-Nitrophenol              | 71. 2,4-Dinitrophenol           | 106. α-Chlordane                    |
| 2. N-Nitrosodimethylamine       | 37. Bis(2-chloroethoxy)methane | 72. 2-Naphthylamine             | 107. Benzidine                      |
| 3. 2-Picoline                   | 38. 2,4-Dichlorophenol         | 73. 2,4-Dinitrotoluene          | 108. Endosulfan II                  |
| 4. N-Nitrosomethylethylamine    | 39. 1,2,4-Trichlorobenzene     | 74. 4-Chlorophenyl phenyl ether | 109. 4,4'-DDE                       |
| 5. 2-Fluorobiphenyl             | 40. Naphthalene-D8             | 75. Fluorene                    | 110. Pyrene                         |
| 6. N-Nitrosodiethylamine        | 41. α-Terpineol                | 76. N-Nitrosodiphenylamine      | 111. p-Terphenyl-D14                |
| 7. Phenol-D5                    | 42. Hexachloro-1,3-butadiene   | 77. Azobenzene                  | 112. Dieldrin                       |
| 8. Phenol                       | 43. Hexachloropropene          | 78. Diphenylamine               | 113. p-Dimethylaminoazobenzene      |
| 9. n-Decane                     | 44. Naphthalene                | 79. 4-Nitroaniline              | 114. Endrin                         |
| 10. Aniline                     | 45. 2,6-Dichlorophenol         | 80. 5-Nitro-o-toluidine         | 115. 4,4'-DDD                       |
| 11. Pentachloroethane           | 46. p-Chloroaniline            | 81. Azobenzene                  | 116. Butyl benzyl phthalate         |
| 12. 2-Chlorophenol              | 47. N-Nitrosodi-n-butylamine   | 82. 2,4,6-Tribromophenol        | 117. Endosulfan II                  |
| 13. Bis(2-chloroethyl)ether     | 48. 4-Chloro-3-methylphenol    | 83. 2-Methyl-4,6-dinitrophenol  | 118. Endrin aldehyde                |
| 14. 1,3-Dichlorobenzene         | 49. 2-Methylnaphthalene        | 84. 4-Bromophenyl phenyl ether  | 119. 3,3'-Dimethylbenzidine         |
| 15. 1,4-Dichlorobenzene-D4      | 50. Hexachlorocyclopentadiene  | 85. n-Octadecane                | 120. 4,4'-DDT                       |
| 16. 1,4-Dichlorobenzene         | 51. 1,2,4,5-Tetrachlorobenzene | 86. α-BHC                       | 121. Endosulfan sulfate             |
| 17. Benzyl alcohol              | 52. 2,3-Dichloroaniline        | 87. Hexachlorobenzene           | 122. 2-Acetylaminofluorene          |
| 18. 1,2-Dichlorobenzene         | 53. 2,4,6-Trichlorophenol      | 88. 4-Aminobiphenyl             | 123. Bis(2-ethylhexyl) phthalate    |
| 19. o-Cresol                    | 54. 2-Fluorophenol             | 89. Pentachlorophenol           | 124. Methoxychlor                   |
| 20. Bis(2-chloroisopropyl)ether | 55. 2,4,5-Trichlorophenol      | 90. γ-BHC                       | 125. Endrin ketone                  |
| 21. p-Cresol                    | 56. 2-Chloronaphthalene        | 91. Pentachloronitrobenzene     | 126. 3,3'-Dichlorobenzidine         |
| 22. m-Cresol                    | 57. 2-Nitroaniline             | 92. Phenanthrene-D10            | 127. Benz[a]anthracene              |
| 23. 3-Methylphenol              | 58. Dimethyl phthalate         | 93. Phenanthrene                | 128. Chrysene-D12                   |
| 24. Hexachloroethane            | 59. Acenaphthylene             | 94. Anthracene                  | 129. Chrysene                       |
| 25. Nitrosodi-n-propylamine     | 60. 2,6-Dinitrotoluene         | 95. Dinoseb                     | 130. Di-n-octyl phthalate           |
| 26. N-Nitroso-di-n-propylamine  | 61. Acenaphthene-D10           | 96. 1,3,5-Trinitrobenzene       | 131. 7,12-Dimethylbenz(a)anthracene |
| 27. o-Toluidine                 | 62. Acenaphthene               | 97. β-BHC                       | 132. Benzo[b]fluoranthene           |
| 28. Acetophenone                | 63. 3-Nitroaniline             | 98. Carbazole                   | 133. Benzo[k]fluoranthene           |
| 29. N-Nitrosomorpholine         | 64. 1,3-Dinitrobenzene         | 99. δ-BHC                       | 134. Hexachlorophene                |
| 30. N-Nitrosopyrrolidine        | 65. 4-Nitrophenol              | 100. Heptachlor                 | 135. Benzo[a]pyrene                 |
| 31. Nitrobenzene-D5             | 66. Pentachlorobenzene         | 101. Di-n-butyl phthalate       | 136. Perylene-D12                   |
| 32. Nitrobenzene                | 67. 1-Naphthylamine            | 102. Aldrin                     | 137. Dibenz[a,h]anthracene          |
| 33. N-Nitrosopiperidine         | 68. Diethyl phthalate          | 103. Heptachlor epoxide         | 138. Indeno[1,2,3-cd]pyrene         |
| 34. Isophorone                  | 69. 2,3,4,6-Tetrachlorophenol  | 104. Fluoranthene               | 139. N-Nitrosodiphenylamine         |
| 35. 2,4-Dimethylphenol          | 70. Diethyl phthalate          | 105. γ-Chlordane                | 140. Benzo[g,h,i]perylene           |




# Environmental Applications


## EPA Method 8081 - Chlorinated Pesticides



## EPA Method 8151 - Chlorinated Herbicides




## Fast H-53 Analysis of Diesel Fuel



**Column:** Zebron ZB-5HT Inferno™  
**Dimensions:** 15 meter x 0.32 mm x 0.10 µm  
**Part No.:** 7EM-G015-02  
**Injection:** On-Column @ 63 °C, 0.1 µL  
**Carrier Gas:** Helium @ 2.7 mL/min (constant flow)  
**Oven Program:** 60 °C to 375 °C @ 25 °C/min  
**Detector:** FID @ 400 °C  
 Note: Diesel Fuel was 200 ppm in dichloromethane with Internal Standards at 50 ppm  
**Sample:** 1. Decane (C10)  
 2. Tetracontane (C40)


## Mineral Oil using H-53 Conditions



**Column:** Zebron ZB-5HT Inferno™  
**Dimensions:** 30 meter x 0.25 mm x 0.10 µm  
**Part No.:** 7HG-G015-02  
**Injection:** On-Column @ 53 °C, 0.1 µL  
**Carrier Gas:** Helium @ 1.3 mL/min (constant flow)  
**Oven Program:** 50 °C for 6 min to 400 °C @ 20 °C/min for 15 min  
**Detector:** FID @ 415 °C  
 Note: Fuel was 10 mg/mL in dichloromethane with 50 ppm markers  
**Sample:** 1. Decane (C10)  
 2. Tetracontane (C40)




## Polynuclear Aromatic Hydrocarbons (PAHs)



**Column:** Zebron ZB-35HT Inferno™  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G025-11

Conditions for both columns:  
**Injection:** Splitless @ 280 °C, 1 µL  
**Carrier Gas:** Helium @ 1.2 mL/min (constant flow)  
**Oven Program:** 80 °C for 0.66 min to 250 °C @ 20 °C/min to 300 °C @ 8 °C/min to 360 °C @ 20 °C/min for 6 min  
**Detector:** MSD @360 °C; 45-400 amu  
 Note: Pulsed splitless injection @ 20 psi for 0.66 min. All analytes @ 10 ppm in dichloromethane.

## Polynuclear Aromatic Hydrocarbons (PAHs)




**Sample:**

|                        |                            |
|------------------------|----------------------------|
| 1. Naphthalene         | 10. Benz[a]anthracene      |
| 2. 2-Methylnaphthalene | 11. Chrysene               |
| 3. Acenaphthylene      | 12. Benzo[b]fluoranthene   |
| 4. Acenaphthene        | 13. Benzo[k]fluoranthene   |
| 5. Fluorene            | 14. Benzo[a]pyrene         |
| 6. Phenanthrene        | 15. Indeno[1,2,3-cd]pyrene |
| 7. Anthracene          | 16. Dibenzo[a,h]anthracene |
| 8. Fluoranthene        | 17. Benzo[g,h,i]perylene   |
| 9. Pyrene              |                            |

# Environmental Applications


## PolyBrominated Diphenyl Ethers (PBDE)

In a GC column, late eluting PBDEs often break down and become difficult to detect. The Zebron™ Inferno™ columns have increased sensitivity for late eluting compounds and thus can easily detect these analytes.



**Column:** Standard 5-type phase  
**Dimensions:** 15 meter x 0.25 mm x 0.10 µm  
**Injection:** On-Column @ 73 °C, 0.5 µL  
**Carrier Gas:** Helium @ 3.4 mL/min (constant pressure)  
**Oven Program:** 70 °C to 160 °C @ 25 °C/min to 320 °C @ 10 °C/min hold 10 min  
**Detector:** ECD @ 340 °C  
**Note:** PBDE Standard @ 2.5 ppm in isooctane


**Sample:** 1. BDE-25 11. BDE-99  
2. BDE-28 12. BDE-85  
3. BDE-75 13. BDE-154  
4. BDE-49 14. BDE-153  
5. BDE-71 15. BDE-138  
6. BDE-47 16. BDE-183  
7. BDE-66 17. BDE-190  
8. BDE-77 18. BDE-203  
9. BDE-100  
10. BDE-119



**Column:** Zebron ZB-5HT Inferno  
**Dimensions:** 15 meter x 0.25 mm x 0.10 µm  
**Part No.:** 7EG-G015-02  
**Injection:** On-Column @ 73 °C, 0.5 µL  
**Carrier Gas:** Helium @ 1.5 mL/min (constant flow)  
**Oven Program:** 70 °C to 160 °C @ 25 °C/min to 350 °C @ 10 °C/min for 5 min  
**Detector:** ECD @ 400 °C  
**Note:** Sample was 2.5 ppm in isooctane.

**Sample:** 1. BDE-25 11. BDE-99  
2. BDE-28 12. BDE-85  
3. BDE-75 13. BDE-154  
4. BDE-49 14. BDE-153  
5. BDE-71 15. BDE-138  
6. BDE-47 16. BDE-183  
7. BDE-66 17. BDE-190  
8. BDE-77 18. BDE-203  
9. BDE-100 19. BDE-209  
10. BDE-119


## EPA Methods 508 and 8081A



**Column:** Zebron ZB-35HT Inferno  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G025-11  
**Injection:** Split 10:1 @ 200 °C, 1 µL  
**Carrier Gas:** Helium @ 0.5 mL/min (constant flow)  
**Oven Program:** 100 °C for 0.5 min to 210 °C @ 30 °C/min for 0.5 min to 230 °C @ 10 °C/min for 2 min to 260 °C @ 5 °C/min for 2 min to 300 °C @ 6 °C/min for 2 min  
**Detector:** MSD @ 350 °C; 50-550 amu  
**Note:** All analytes were 10 ppm in isooctane.

|                              |                             |                        |
|------------------------------|-----------------------------|------------------------|
| 1. DBCP                      | 15. δ-BHC                   | 29. Nitrofen           |
| 2. Hexachlorocyclopentadiene | 16. Chlorothalonil          | 30. Kepone             |
| 3. Etridiazole               | 17. Aldrin                  | 31. 4,4'-DDD           |
| 4. Chloroneb                 | 18. DCPA (Dacthal)          | 32. Endosulfan II      |
| 5. Propachlor                | 19. Heptachlor epoxide      | 33. 4,4'-DDT           |
| 6. Diallyl                   | 20. trans-Chlordane (gamma) | 34. Endrin Aldehyde    |
| 7. Hexachlorobenzene         | 21. trans-Nonachlor         | 35. Endosulfan sulfate |
| 8. α-BHC                     | 22. cis-Chlordane (alpha)   | 36. Captafol           |
| 9. PCNB                      | 23. Endosulfan I            | 37. Methoxychlor       |
| 10. γ-BHC                    | 24. 4,4'-DDE                | 38. Endrin Ketone      |
| 11. β-BHC                    | 25. Dieldrin                | 39. Permethrin         |
| 12. Alachlor                 | 26. Chloropropylate         | 40. Permethrin isomer  |
| 13. Heptachlor               | 27. Chlorobenzilate         |                        |
| 14. Dichlone                 | 28. Endrin                  |                        |


## Poly Brominated Diphenyl Ethers (PBDEs)




**Column:** Zebron ZB-5HT Inferno  
**Dimensions:** 30 meter x 0.25 mm x 0.10 µm  
**Part No.:** 7HG-G015-02  
**Injection:** Pulsed @ 350 °C, 2 µL  
**Carrier Gas:** Helium @ 1.3 mL/min (constant flow)  
**Oven Program:** 125 °C to 380 °C @ 22 °C/min for 15 min  
**Detector:** MSD; 45 - 800 amu  
**Note:** Analytes were 2.5 ppm in isooctane. Pulsed splitless injection at 50 psi using direct connect liner

|           |             |             |
|-----------|-------------|-------------|
| 1. BDE-25 | 8. BDE-77   | 15. BDE-138 |
| 2. BDE-28 | 9. BDE-100  | 16. BDE-183 |
| 3. BDE-75 | 10. BDE-119 | 17. BDE-190 |
| 4. BDE-49 | 11. BDE-99  | 18. BDE-203 |
| 5. BDE-71 | 12. BDE-85  |             |
| 6. BDE-47 | 13. BDE-154 |             |
| 7. BDE-66 | 14. BDE-153 |             |

## Motor Oil



## Paraffin Wax




## ASTM Methods D2887 and D6352 on One Column

ASTM D2887 and D6352 both set standard test methods for boiling range determination of petroleum fractions by GC. While ASTM D2887 is used for products having a final boiling point of 538 °C or lower, ASTM D6352 is used for determining petroleum distillate fractions with a boiling point (BP) up to 700 °C. The Zebron™ Inferno™ columns make

simulated distillation analysis easier and less expensive. Typically you would have to use a traditional fused silica column for ASTM D2887 and a metal column for ASTM D6352. Zebron Inferno columns provide a robust solution for both analyses.

## ASTM Method D6352




**Column:** Zebron ZB-1HT Inferno™  
**Dimensions:** 5 meter x 0.53 mm x 0.10 µm  
**Part No.:** 7AK-G014-02  
**Injection:** On-Column @ 43 °C, 0.1 µL  
**Carrier Gas:** Helium @ 4.4 mL/min (constant flow)  
**Oven Program:** 40 °C for 0.5 min to 430 °C @ 20 °C/min for 10 min  
**Detector:** FID @ 430 °C

**Note:** Sample was a combination of PolyWax 655 and retention time markers C8-C40 in CS2/Chloroform.

|             |             |         |         |         |         |
|-------------|-------------|---------|---------|---------|---------|
| 1. C10      | 11. Phytane | 21. C28 | 31. C38 | 41. C56 | 51. C76 |
| 2. C11      | 12. C19     | 22. C29 | 32. C39 | 42. C58 | 52. C78 |
| 3. C12      | 13. C20     | 23. C30 | 33. C40 | 43. C60 | 53. C80 |
| 4. C13      | 14. C21     | 24. C31 | 34. C42 | 44. C62 | 54. C82 |
| 5. C14      | 15. C22     | 25. C32 | 35. C44 | 45. C64 | 55. C84 |
| 6. C15      | 16. C23     | 26. C33 | 36. C46 | 46. C66 | 56. C86 |
| 7. C16      | 17. C24     | 27. C34 | 37. C48 | 47. C68 | 57. C88 |
| 8. C17      | 18. C25     | 28. C35 | 38. C50 | 48. C70 | 58. C90 |
| 9. Pristane | 19. C26     | 29. C36 | 39. C52 | 49. C72 |         |
| 10. C18     | 20. C27     | 30. C37 | 40. C54 | 50. C74 |         |

## ASTM Method D2887-97



**Column:** Zebron ZB-1HT Inferno™  
**Dimensions:** 15 meter x 0.53 mm x 0.15 µm  
**Part No.:** 7EK-G014-05  
**Injection:** On-Column @ 38 °C, 0.1 µL  
**Carrier Gas:** Helium @ 10 mL/min (constant flow)  
**Oven Program:** 35 °C for 1 min to 360 °C @ 25 °C/min  
**Detector:** FID @ 375 °C


|                |                       |
|----------------|-----------------------|
| 1. Hexane      | 10. Octadecane        |
| 2. Heptane     | 11. Eicosane          |
| 3. Octane      | 12. Tetracosane       |
| 4. Nonane      | 13. Octacosane        |
| 5. Decane      | 14. n-Dotriacontane   |
| 6. Undecane    | 15. n-Hexatriacontane |
| 7. Dodecane    | 16. Tetracontane      |
| 8. Tetradecane | 17. Tetraoctacontane  |
| 9. Hexadecane  |                       |



## Free and Total Glycerin in Biodiesel


The Zebron™ ZB-1HT Inferno™ column is ideal for biodiesel analysis because it is stable at the elevated temperatures needed to resolve and remove high boiling triglycerides. While standard polyimide resin degrades at temperatures above 380 °C, the Inferno columns utilize a polyimide resin material that shows minimal thermal breakdown, even at programmed temperatures up to 430 °C– the highest thermal stability of any non-metal GC column!

## ASTM Method D6584: Free and Total Glycerin in B-100 Biodiesel Methyl Esters




**Column:** Zebron ZB-5HT Inferno™  
**Dimensions:** 15 meter x 0.32 mm x 0.10 µm + 2 meter x 0.53 mm Z-Guard  
**Part No.:** 7EM-G015-02  
**Injection:** On-Column @ 53 °C, 1 µL  
**Carrier Gas:** Helium @ 3.0 mL/min (constant flow)  
**Oven Program:** 50 °C for 1 min to 180 °C @ 15 °C/min to 230 °C @ 7 °C/min to 380 °C @ 30 °C/min for 10 min  
**Detector:** FID @ 380 °C  
 Note: A 2.0 m x 0.53 mm Guard Column was connected to the analytical column per ASTM method requirements.  
**Sample:**

1. Glycerol
2. Butanetriol (ISTD1)
3. Esters
4. Monoglycerides
5. 1-Monooleoyl-rac-glycerol
6. Tricarpin (ISTD2)
7. Diglycerides
8. 1,3-Diolein
9. Triglycerides
10. Triolein



**Column:** Zebron ZB-35HT Inferno™  
**Dimensions:** 15 meter x 0.25 mm x 0.10 µm  
**Part No.:** 7EG-G025-02  
**Injection:** On-Column @ 53 °C, 0.1 µL  
**Carrier Gas:** Helium @ 1.3 mL/min (constant flow)  
**Oven Program:** 50 °C for 1 min to 180 °C @ 15 °C/min to 230 °C @ 7 °C/min to 380 °C @ 30 °C/min for 10 min  
**Detector:** FID @ 380 °C  
 Note: A 30 cm x 0.5 mm ID guard column was connected to the analytical column per ASTM method requirements  
**Sample:**

1. Butanetriol (ISTD1)
2. Ester
3. Monoglycerides
4. 1-Monooleoyl-rac-glycerol
5. Tricarpin (ISTD2)
6. Diglycerides
7. 1,3-Diolein
8. Triglycerides
9. Triolein




**Column:** Zebron ZB-XLB-HT Inferno™  
**Dimensions:** 15 meter x 0.25 mm x 0.10 µm  
**Part No.:** 7HG-G024-11  
**Injection:** On-Column @ 53 °C, 0.10 µL  
**Carrier Gas:** Helium @ 1.3 mL/min (constant flow)  
**Oven Program:** 50 °C for 1 min to 180 °C @ 15 °C/min to 230 °C @ 7 °C/min to 380 °C @ 30 °C/min for 10 min  
**Detector:** FID @ 380 °C  
 Note: A 30 cm x 0.5 mm ID guard column was connected to the analytical column per ASTM method requirements  
**Sample:**

1. Butanetriol (ISTD1)
2. Ester
3. Monoglycerides
4. 1-Monooleoyl-rac-glycerol
5. Tricarpin (ISTD2)
6. Diglycerides
7. 1,3-Diolein
8. Triglycerides
9. Triolein




## Melamine in Baby Formula




**Column:** Zebron ZB-XLB-HT Inferno™  
**Dimensions:** 15 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7EG-G024-11  
**Injection:** On-Column @ 103 °C, 1 µL  
**Carrier Gas:** Helium @ 1.4 mL/min (constant flow)  
**Oven Program:** 100 °C for 0.5 min to 320 °C @ 25 °C/min  
**Detector:** Mass Selective (MSD) @ 325 °C  
**Sample:** Analytes are 200 ng / 100 µL in BSTFA / Pyridine (1:1)  
 1. Cyanuric Acid 13C3 (IS)  
 2. Cyanuric Acid  
 3. Melamine 13C3 15N3 (IS)  
 4. Melamine

## Butter Triglycerides




**Column:** Zebron ZB-5HT Inferno™  
**Dimensions:** 15 meter x 0.32 mm x 0.10 µm  
**Part No.:** 7EM-G015-02  
**Injection:** On-Column @ 103 °C, 2 µL  
**Carrier Gas:** Helium @ 1.8 mL/min (constant flow)  
**Oven Program:** 100 °C to 400 °C @ 14 °C/min for 10 min  
**Detector:** FID @ 400 °C  
**Sample:** Butter

## Sterols from Olive Oil: Underivatized



**Column:** Zebron ZB-5HT Inferno™  
**Dimensions:** 30 meter x 0.25 mm x 0.10 µm  
**Part No.:** 7HG-G015-02  
**Injection:** Splitless @ 350 °C, 0.5 µL  
**Carrier Gas:** Helium @ 2 mL/min (constant flow)  
**Oven Program:** 220 °C to 350 °C @ 15 °C/min  
**Detector:** FID @ 350 °C  
**Sample:**  
 1. Campesterol  
 2. Stigmasterol  
 3. β-Sitosterol  
 4. Betulin (IS)


## Sterols from Olive Oil: Derivatized



**Column:** Zebron ZB-5HT Inferno™  
**Dimensions:** 30 meter x 0.25 mm x 0.10 µm  
**Part No.:** 7HG-G015-02  
**Injection:** Splitless @ 350 °C, 0.5 µL  
**Carrier Gas:** Helium @ 2 mL/min (constant flow)  
**Oven Program:** 220 °C to 350 °C @ 15 °C/min.  
**Detector:** FID @ 350 °C  
**Note:** Sterols derivatized via BSTFA:TMCS; 99:1 in pyridine  
**Sample:**  
 1. Cholesterol  
 2. Campesterol  
 3. Stigmasterol  
 4. β-Sitosterol  
 5. Betulin (IS)

# Chemical Applications

## Amine and Nitrosamine




App ID 18773

**Column:** Zebron ZB-35HT Inferno™  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G025-11  
**Injection:** Split 50:1 @ 350 °C, 1 µL  
**Carrier Gas:** Helium @ 2.1 mL/min (constant flow)  
**Oven Program:** 140 °C to 400 °C @ 10 °C/min  
**Detector:** FID @ 410 °C

Note: Application courtesy of Northeastern Chemical Company

**Sample:** 1. MW 135 Amine  
 2. MW 163 Amine  
 3. MW 178 Amine  
 4. MW 208 Amine  
 5. MW 164 Nitrosamine  
 6. MW 178 Nitrosamine  
 7. MW 192 Nitrosamine  
 8. MW 208 Nitrosamine

## Diethylene Glycol Divinyl Ether



App ID 18770


**Column:** Zebron ZB-35HT Inferno™  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G025-11  
**Injection:** Split 40:1 @ 250 °C, 1 µL  
**Carrier Gas:** Helium @ 1.4 mL/min (constant flow)  
**Oven Program:** 40 °C for 0.5 min to 150 °C @ 10 °C/min  
 to 310 °C @ 25 °C/min for 2 min  
**Detector:** FID @ 320 °C

Note: Application courtesy of Northeastern Chemical Company

**Sample:** 1. Isopar C  
 2. Diethylene glycol divinyl ether  
 3. Isobornyl acrylate

# Pharmaceutical/Drugs

## Drug Screen



App ID 16227

**Column:** Zebron ZB-1HT Inferno™  
**Dimensions:** 30 meter x 0.25 mm x 0.10 µm  
**Part No.:** 7HG-G014-02  
**Injection:** Split 15:1 @ 250 °C, 1 µL  
**Carrier Gas:** Helium 1.1 mL/min (constant flow)  
**Oven Program:** 50 °C for 1 min to 320 °C @ 25 °C/min for 3 min  
**Detector:** MSD @ 310 °C; 45-450 amu

Note: Analytes were 25 ppm

**Sample:** 1. Nicotine  
 2. Contaminant  
 3. Acetaminophen  
 4. Allobarbitol  
 5. Aprobital  
 6. Ibuprofen  
 7. Butabarbital  
 8. Phenacetin  
 9. Amobarbital  
 10. Pentobarbital  
 11. Meprobamate  
 12. Secobarbital  
 13. Caffeine  
 14. Methyl Benzilate  
 15. Benzphetamine  
 16. Hexobarbital  
 17. Dimenhydrinate  
 18. Doxylamine  
 19. Phenobarbital  
 20. Methapyrilene  
 21. Chlorpheniramine  
 22. Brompheniramine  
 23. Cocaine  
 24. Chlorcyclizine  
 25. Codeine  
 26. Diazepam  
 27. Hydrocodone  
 28. Oxymorphone

## Cool-Lock™ Nut

### Capillary GC Installation with the Cool-Lock™ Nut

- Avoid burning your fingers again – cools with the oven
- No Wite-Out® or Tipp-Ex® needed
- Achieve the proper installation depth each and every time
- No need for wrench with hand-tightened connections



- 1** Insert column through column nut and ferrule. Insert the end of the column and column nut into the bottom of the Cool-Lock Nut Installation Gauge.



- 2** Hand-tighten the column nut into the gauge until the ferrule starts to compress. Adjust the length of the column to the desired length.



- 3** Unscrew the installation gauge and install the nut and column into the instrument.

- 4** Make sure Cool-Lock nut is in the full extension position to achieve the correct pre-measured position, and allow for proper leak checking.



#### Ordering Information

##### Cool-Lock GC Capillary Nut For Agilent GC Systems\*

| Part No. | Description                                                  | Unit | Price |
|----------|--------------------------------------------------------------|------|-------|
| AGO-8319 | Cool-Lock GC Capillary Nut For Use With Short-Style Ferrules | ea   |       |
| AGO-8320 | Cool-Lock GC Capillary Nut For Use With Long-Style Ferrules  | ea   |       |
| AGO-8349 | Cool-Lock Nut Installation Gauge                             | ea   |       |

#### Replacement Ferrules

|          |                                                      |       |  |
|----------|------------------------------------------------------|-------|--|
| AGO-4701 | GC Capillary Ferrules Graphite 1/16 in. to 0.5 mm ID | 10/pk |  |
| AGO-4704 | GC Capillary Ferrules Graphite 1/16 in. to 0.8 mm ID | 10/pk |  |

\*Guaranteed fit for Agilent 5850, 5890, 6850, 6890 GC systems



#### Cool-Lock GC Capillary Nut For Shimadzu Systems\*\*

| Part No. | Description                                                  | Unit | Price |
|----------|--------------------------------------------------------------|------|-------|
| AGO-8419 | Cool-Lock GC Capillary Nut For Use With Short-Style Ferrules | ea   |       |
| AGO-8420 | Cool-Lock Nut Installation Gauge                             | ea   |       |

\*\*Guaranteed fit for Shimadzu 2010 and 2014 GC systems



## Guardian™: Integrated Guard Columns

- Eliminates the potential for leaks
- Extends column life
- Improves analyte focusing for low boiling compounds
- Aggressively tested to ensure deactivation

### Traditional Union Connected Guard Column

vs.

### Guardian Integrated Guard Column



#### Ordering Information

##### Guardian: Integrated Guard Column

| GC Phase      | Dimensions                        | 5 m Guardian Part No. | Price |
|---------------|-----------------------------------|-----------------------|-------|
| Zebron ZB-1HT | 30 meter x 0.25 mm x 0.10 df (µm) | 7HG-G014-02-GGA       |       |
| Zebron ZB-5HT | 30 meter x 0.25 mm x 0.10 df (µm) | 7HG-G015-02-GGA       |       |
| Zebron ZB-5HT | 30 meter x 0.25 mm x 0.25 df (µm) | 7HG-G015-11-GGA       |       |

## Inlet Base Seals

- Prevents analyte adsorption
- Improves detection sensitivity
- Improves method reproducibility
- Fits all Agilent 5890/6890 split/splitless injection ports



### Ordering Information

#### Gold Inlet Base Seals

| Part No.                                                                            | Description                                                    | Similar to Mfr. No.* | Unit  | Price |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------|----------------------|-------|-------|
| <b>Standard, single groove for splitless applications, 0.8 mm dia. inlet hole</b>   |                                                                |                      |       |       |
| AGO-7518                                                                            | Gold Inlet Base Seal, splitless (single groove), 0.8 mm        | 18740-20885          | 2/pk  |       |
| AGO-7519                                                                            | Gold Inlet Base Seal, splitless (single groove), 0.8 mm        | 18740-20885          | 10/pk |       |
| <b>High Split Flow, cross groove for split applications, 0.8 mm dia. inlet hole</b> |                                                                |                      |       |       |
| AGO-7520                                                                            | Gold Inlet Base Seal, split (double groove/cross), 0.8 mm      | 5182-9652            | 2/pk  |       |
| AGO-7521                                                                            | Gold Inlet Base Seal, split (double groove/cross), 0.8 mm      | 5182-9652            | 10/pk |       |
| <b>Replacement Inlet Seal Washers</b>                                               |                                                                |                      |       |       |
| AGO-8397                                                                            | Inlet Seal Washers, Gold plated, for Agilent GC injection port | -                    | 12/pk |       |



### Stainless Steel Base Seals

| Part No.                                                                            | Description                                                          | Similar to Mfr. No.* | Unit  | Price |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------|----------------------|-------|-------|
| <b>Standard, single groove for splitless applications, 0.8 mm dia. inlet hole</b>   |                                                                      |                      |       |       |
| AGO-8393                                                                            | Stainless Steel Inlet Base Seal, splitless (single groove), 0.8 mm   | 18740-20880          | 2/pk  |       |
| AGO-8394                                                                            | Stainless Steel Inlet Base Seal, splitless (single groove), 0.8 mm   | 18740-20880          | 10/pk |       |
| <b>High Split Flow, cross groove for split applications, 0.8 mm dia. inlet hole</b> |                                                                      |                      |       |       |
| AGO-8395                                                                            | Stainless Steel Inlet Base Seal, split (double groove/cross), 0.8 mm | -                    | 2/pk  |       |
| AGO-8396                                                                            | Stainless Steel Inlet Base Seal, split (double groove/cross), 0.8 mm | -                    | 10/pk |       |
| <b>Replacement Inlet Seal Washers</b>                                               |                                                                      |                      |       |       |
| AGO-7522                                                                            | Inlet Seal Washers, Stainless Steel, for Agilent GC injection port   | 5061-5869            | 12/pk |       |

## Injection Port Inlet Liners

- Fit standard instrument manufacturers' inlets
- Popular designs at excellent prices
- Reduce analyte discrimination
- Improve reproducibility and results
- Clean, non-adsorptive/non-reactive surface



### Ordering Information

| Description                                         | GC Model No.   | Dimensions ID x L x OD (mm) | Material* (deactivated) | Quartz Wool (Y/N) | Similar to Manufacturer Number** | Part No. | Unit | Price |
|-----------------------------------------------------|----------------|-----------------------------|-------------------------|-------------------|----------------------------------|----------|------|-------|
| <b>Agilent Technologies (HP)</b>                    |                |                             |                         |                   |                                  |          |      |       |
| Single Taper Direct Connect with Side Hole (top)    | 5880/5980/6890 | 4 x 78.5 x 6.3              | B (y)                   | N                 | 21055                            | AGO-7850 | 5/pk |       |
| Single Taper Direct Connect with Side Hole (bottom) | 5880/5980/6890 | 4 x 78.5 x 6.3              | B (y)                   | N                 | G1544-80730 20771                | AGO-7851 | 5/pk |       |
| Split / Splitless Focus Liner with wool             | 5880/5890/6890 | 4 x 78.5 x 6.3              | B (y)                   | Y                 | 210-4004                         | AGO-7515 | 5/pk |       |
| <b>Shimadzu</b>                                     |                |                             |                         |                   |                                  |          |      |       |
| Splitless Liner                                     | 9A/15A/16A     | 3.4 x 139 x 5.0             | B (y)                   | N                 | 20749                            | AGO-4669 | 5/pk |       |
| <b>Thermo Scientific (Finnigan)</b>                 |                |                             |                         |                   |                                  |          |      |       |
| Splitless / Single Taper Liner /Trace               | 8000 series    | 5 x 105 x 8.0               | B (y)                   | N                 | 45350033                         | AGO-7852 | 5/pk |       |

\* B= Borosilicate; Deactivated = Yes (y) or No (n). \*\* Similar to but not always an exact equivalent to the original manufacturer's product.



Additional GC accessories are available.  
For more information contact our GC Specialists at [GCspecialist@phenomenex.com](mailto:GCspecialist@phenomenex.com).

## Septa

### Silicone Rubber Inlet Septa

- Advanced silicone formulations reduce coring and increase septum lifetime
- Enhanced durability and re-sealing capabilities
- Precision formed and cut to exacting specifications
- Reduced needle damage
- General purpose and specialty low bleed, high-temperature septa available

### PhenoRed™-400

Fire brick red, with a 50 durometer rating, this specially formulated high-temperature, low-bleed silicone rubber septum is designed and conditioned for high sensitivity analysis where the demands for low bleed and heat stability are critical. Say goodbye to ghost peaks! Temperature-rated to 400 °C.

### PhenoGreen™-400

Long life, high-temperature green septum, rated to 400 °C. GuideRight™ through-hole to guide the needle during injection.

## GuideRight™ Injection Hole Septa

The GuideRight™ through-hole guides the needle during injection. Septum performance and lifetime increase while downtime due to bent needles or dirty injection ports decreases. Choose from either high-temperature PhenoGreen™ or PhenoRed™, both rated to 400 °C.

## Long Ferrules

### Graphite Ferrules

- High-purity graphite significantly reduces ferrule bleed
- Special construction minimizes "flaking"
- Stable to 450 °C - excellent for high temperature use
- Not for use with GC/MS transfer lines



#### Ordering Information Silicone Rubber Septa

| Part No.                                   | Description                 | Unit   | Price |
|--------------------------------------------|-----------------------------|--------|-------|
| <b>PhenoRed 3/16 in. (9.5 mm) Diameter</b> |                             |        |       |
| AGO-4690                                   | PhenoRed-400 Injector Septa | 50/pk  |       |
| AGO-4691                                   | PhenoRed-400 Injector Septa | 100/pk |       |
| <b>PhenoRed 7/16 in. (11 mm)</b>           |                             |        |       |
| AGO-4696                                   | PhenoRed-400 Injector Septa | 50/pk  |       |
| AGO-4697                                   | PhenoRed-400 Injector Septa | 100/pk |       |

#### Ordering Information GuideRight Injection Hole Septa

| Part No.                                            | Description                    | Unit  | Price |
|-----------------------------------------------------|--------------------------------|-------|-------|
| <b>PhenoGreen 3/16 in. (9.5 mm) Diameter 400 °C</b> |                                |       |       |
| AGO-7874                                            | PhenoGreen™ 400 Injector Septa | 50/pk |       |
| <b>PhenoGreen 7/16 in. (11 mm) Diameter 400 °C</b>  |                                |       |       |
| AGO-7875                                            | PhenoGreen™ 400 Injector Septa | 50/pk |       |
| <b>PhenoRed 3/16 in. (9.5 mm) Diameter 400 °C</b>   |                                |       |       |
| AGO-7916                                            | PhenoRed™ 400 Injector Septa   | 50/pk |       |
| <b>PhenoRed 7/16 in. (11 mm) Diameter 400 °C</b>    |                                |       |       |
| AGO-7917                                            | PhenoRed™ 400 Injector Septa   | 50/pk |       |

#### Ordering Information Graphite Ferrules, Rated to 450 °C

| Part No.                 | Description                         | Unit  | Price |
|--------------------------|-------------------------------------|-------|-------|
| <b>0.4 mm Ferrule ID</b> |                                     |       |       |
| AGO-4698                 | Graphite Ferrule 1/16 in. to 0.4 mm | 10/pk |       |
| AGO-4699                 | Graphite Ferrule 1/16 in. to 0.4 mm | 50/pk |       |
| <b>0.5 mm Ferrule ID</b> |                                     |       |       |
| AGO-4701                 | Graphite Ferrule 1/16 in. to 0.5 mm | 10/pk |       |
| AGO-4702                 | Graphite Ferrule 1/16 in. to 0.5 mm | 50/pk |       |
| <b>0.8 mm Ferrule ID</b> |                                     |       |       |
| AGO-4704                 | Graphite Ferrule 1/16 in. to 0.8 mm | 10/pk |       |
| AGO-4705                 | Graphite Ferrule 1/16 in. to 0.8 mm | 50/pk |       |



Additional GC accessories are available.  
For more information contact our GC Specialists at [GCspecialist@phenomenex.com](mailto:GCspecialist@phenomenex.com).

#### Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

#### Trademarks

Zebron, Inferno, Cool-Lock, Guardian, Engineered Self Cross-linking, ESC, PhenoGreen, PhenoRed, and GuideRight are trademarks of Phenomenex, Inc. Cool-Lock Nut is patent pending by Phenomenex. Agilent and J&W are registered trademarks of Agilent Technologies, Inc. Varian is a registered trademark under license to Varian, Inc. Restek and Rtx are registered trademarks of Restek Corporation. Mobil is a registered trademark of Exxon Mobil Corporation. Wite-Out and TIPP-EX are registered trademarks of BIC Corporation.

#### Disclaimer

Comparative separations may not be representative of all applications.

Cool-Lock Nut is U.S. patent pending by Phenomenex, Inc.

© 2010 Phenomenex, Inc. All rights reserved.

guarantee

If Zebron does not provide you with equivalent separations as compared to any other GC column of the same phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.



#### Australia

t: 02-9428-6444  
f: 02-9428-6445  
auinfo@phenomenex.com

#### Austria

t: 01-319-1301  
f: 01-319-1300  
anfrage@phenomenex.com

#### Belgium

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

#### Canada

t: (800) 543-3681  
f: (310) 328-7768  
info@phenomenex.com

#### Denmark

t: 4824 8048  
f: 4810 6265  
nordicinfo@phenomenex.com

#### Finland

t: 09-4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

#### France

t: 01 30 09 21 10  
f: 01 30 09 21 11  
franceinfo@phenomenex.com

#### Germany

t: 06021-58830-0  
f: 06021-58830-11  
anfrage@phenomenex.com

#### Ireland

t: 01 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

#### Italy

t: 051 6327511  
f: 051 6327555  
italiainfo@phenomenex.com

#### Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

#### Mexico

t: (55) 5018 3791  
f: (310) 328-7768  
tecnicomx@phenomenex.com

#### Netherlands

t: 030-2418700  
f: 030-2383749  
nlinfo@phenomenex.com

#### New Zealand

t: 09-4780951  
f: 09-4780952  
nzinfo@phenomenex.com

#### Norway

t: 81 00 20 05  
f: +45 4810 6265  
nordicinfo@phenomenex.com

#### Puerto Rico

t: (800) 541-HPLC  
f: (310) 328-7768  
info@phenomenex.com

#### United Kingdom

t: 01625-501367  
f: 01625-501796  
ukinfo@phenomenex.com

All other countries:  
Corporate Office USA



t: (310) 212-0555  
f: (310) 328-7768  
info@phenomenex.com



[www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)