

La-Pha-Pack®

Werner Reifferscheidt GmbH

Consumables for Waters



La-Pha-Pack® Werner Reifferscheidt GmbH, Am Parir 20, 52379 Langerwehe, Internet: www.la-pha-pack.com
Tel.: (+49) 24 23 / 94 31-0, Fax: (+49) 24 23 / 94 31-33, E-Mail: Service@La-Pha-Pack.com

La-Pha-Pack® • Werner Reifferscheidt GmbH • Tel.: (+49) 24 23 / 94 31-0 • Fax: (+49) 24 23 / 94 31-33 • E-Mail: Service@La-Pha-Pack.com

PP MICRO-VIALS ND9 + ND11



Consumables for Waters

Waters Alliance®
Alliance® 2690/2695
Alliance® 2790/2795
ACQUITY UPLC™
Waters Alliance® HT System

MICROLITER ND9 + ND11



Short Thread ND9									
1:1 Drawings									
LPP Art.-No.	11 09 0500 (11 09 1241)	11 09 0519	11 09 0520 (11 09 1242)	11 19 1205	11 19 1516	11 09 0999	11 09 0620	11 14 1189	
See Waters Art.-No.	186000273 (186000273DV)	186000848	186000848DV	Alt. 186002634		186002804	186002802		
LPP Art.-Description	1.5ml Short Thread Vial, 32 x 11.6mm, clear glass, 1" hydrol. class, wide opening <i>(SILANIZED)</i>	1.5ml Short Thread Vial, 32 x 11.6mm, clear glass, 1" hydrol. class, wide opening, label + filling lines	1.5ml Short Thread Vial, 32 x 11.6mm, amber glass, 1" hydrol. class, wide opening, label + filling lines <i>(SILANIZED)</i>	1.5ml Polypropylene Short Thread Vial, 32 x 11.6mm, transp., with filling lines, (slightly konkave interior bottom shape) <i>(SILANIZED)</i>	1.5ml Polypropylene Short Thread Vial, 32 x 11.6mm, amber, with filling lines, (slightly konkave interior bottom shape)	Short Thread Vial with integrated 0.2ml Micro-Insert, clear glass, 1" hydrol. class, label + filling lines	1.1ml Micro-Insert Short Thread Vial ND9, 32 x 11.6mm, clear glass, 1" hydrol. class	TPX Short Thread Vial ND9, 32 x 11.6mm, with integrated 0.2ml Glass Micro-Insert	Top Serr™ Micro-Insert
LPP Packaging Unit	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PE-Bag	100 pieces per PE-Bag	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box

Crimp Neck ND11					PP Vials ND9 + ND11				
1:1 Drawings									
LPP Art.-No.	11 09 0356	11 09 0476	11 09 0477	11 09 0921	11 09 0619	11 14 1190	11 19 0932	11 19 0933	
See Waters Art.-No.	WAT094222		WAT094223				186002626	186002628	
LPP Art.-Description	1.5ml Crimp Neck Vial, 32 x 11.6mm, clear glass, 1" hydrol. class, wide opening	1.5ml Crimp Neck Vial, 32 x 11.6mm, clear glass, 1" hydrol. class, wide opening, label + filling lines	1.5ml Crimp Neck Vial, 32 x 11.6mm, amber glass, 1" hydrol. class, wide opening, label + filling lines	Crimp Neck Vial with integrated 0.2ml Micro-Insert, clear glass, 1" hydrol. class, label + filling lines	1.1ml Microliter Crimp Neck Vial ND11, 32 x 11.6mm, clear glass, 1" hydrol. class	Top Serr™ TPX Snap Ring Vial ND11, 32 x 11.6mm, with integrated 0.2ml Glass Micro-Insert	0.3ml Polypropylene Short Thread Micro-Vial ND9, 32 x 11.6mm, transparent <i>(also available in TPX or amber PP)</i>	0.3ml Polypropylene Snap Ring Micro-Vial ND9, 32 x 11.6mm, transparent <i>(also available in TPX or amber PP)</i>	
LPP Packaging Unit	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PE-Bag	100 pieces per PE-Bag	

HPLC

2




Consumables for Waters

Waters Alliance®
Alliance® 2690/2695
Alliance® 2790/2795
ACQUITY UPLC™
Waters Alliance® HT System



LPP Art.-No.	11 09 0627	11 09 0644	11 09 0645	11 14 1190	10 09 0743	10 09 1196	10 09 1197
See Waters Art.-No.	WAT094219		WAT094220	WAT063300			
LPP Art.-Description	1.5ml Snap Ring Vial, 32 x 11.6mm, clear glass, 1" hydrof. class, wide opening	1.5ml Snap Ring Vial, 32 x 11.6mm, clear glass, 1" hydrof. class, wide opening, label + filling lines	1.5ml Snap Ring Vial, 32 x 11.6mm, amber glass, 1" hydrof. class, wide opening, label + filling lines	TopSert™ TPX Snap Ring Vial ND11, 32 x 11.6mm, with integrated 0.2ml Glass Micro-Insert	1.5ml Screw Neck Vial ND10, 32 x 11.6mm, clear glass, 1" hydrof. class, wide opening	1.5ml Screw Neck Vial ND10, 32 x 11.6mm, clear glass, 1" hydrof. class, wide opening, label + filling lines	1.5ml Screw Neck Vial ND10, 32 x 11.6mm, amber glass, 1" hydrof. class, wide opening, label + filling lines
LPP Packaging Unit	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box

MICRO-INSERTS FOR WIDE OPENINGS (Short Thread ND9, Crimp Neck ND11, Snap Ring ND11, Screw Neck ND10)



LPP Art.-No.	06 09 0357	06 09 0669	06 09 0866	06 09 0865 (06 09 1343)	06 09 1240
See Waters Art.-No.	WAT094170 (WAT094170DV)				
LPP Art.-Description	0.1ml Micro-Insert, 31 x 6mm, clear glass, 15mm top	0.1ml Micro-Insert, 31 x 6mm, clear glass, 12mm top	0.1ml Micro-Insert, 31 x 6mm, clear glass, flat bottom	0.1ml Micro-Insert, 29 x 5.7mm, clear glass, with assembled plastic spring (SILANIZED)	0.1ml Micro-Insert, 31 x 6mm, clear glass, 15mm top SILANIZED
LPP Packaging Unit	10 x 100 pieces per PP-Box	10 x 100 pieces per PP-Box	10 x 100 pieces per PP-Box	10 x 100 pieces per PP-Box	10 x 100 pieces per PP-Box

Crimpers

11 06 0006	11 07 0002
PSL904301	
11mm Crimper	11mm Decapper
1 piece per PP-Box	1 piece per PP-Box



Vial Racks

12 21 1001
Vial Rack, Acrylic (173 x 95 x 20mm), 50 cavities with a diameter of 12mm
1 piece



HPLC

3



Consumables for Waters

Waters Alliance®
Alliance® 2690/2695
Alliance® 2790/2795
ACQUITY UPLC™
Waters Alliance® HT System



Seals for Short Thread ND9 (with assembled Septum)

09 15 0982		9mm PP Short Thread Cap, blue, centre hole;	PTFE virginal, 53° shore D, 0.2mm
09 15 0867		9mm PP Short Thread Cap, blue, centre hole;	Natural Rubber red-orange/TEF transp., 60° shore A, 1.0mm
09 15 0838	alt. 186000274	9mm PP Short Thread Cap, blue, centre hole;	UltraClean Silicone white/PTFE red, 55° shore A, 1.0mm
09 15 0868		9mm PP Short Thread Cap, blue, centre hole;	PTFE red/Silicone white/PTFE red, 45° shore A, 1.0mm
09 15 0869	alt. 186000305	9mm PP Short Thread Cap, blue, centre hole;	Silicone white/PTFE blue, with slit, 55° shore A, 1.0mm
09 04 1220	alt. 186000274	9mm UltraBond™ Seal: PP Short Thread Cap, black, centre hole;	Silicone white/PTFE red, 45° shore A, 1.3mm
09 04 1533	186000274	9mm UltraBond™ Seal: PP Short Thread Cap, blue, centre hole;	Silicone white/PTFE beige, 45° shore A, 1.3mm
09 04 1534	186000305	9mm UltraBond™ Seal: PP Short Thread Cap, blue, centre hole;	Silicone white/PTFE beige, with slit, 45° shore A, 1.3mm Improved WATERS LectraBond Quality!

Further 9mm Short Thread Seals with different Septa resp. with red, transparent, green, yellow or black caps are available!
Packaging Unit: 100 pieces per PE-Bag

Seals for Crimp Neck ND11 (with assembled Septum)

11 03 0300	PSL412115	11mm Alu. Crimp Cap, clear lacquered, centre hole;	Natural Rubber red-orange/Butyl red/TEF transp., 45° shore A, 1.0mm
11 03 0209		11mm Alu. Crimp Cap, clear lacquered, centre hole;	Natural Rubber red-orange/TEF transparent, 60° shore A, 1.0mm
11 03 0339	PSL408027	11mm Alu. Crimp Cap, clear lacquered, centre hole;	PTFE virginal, 53° shore A, 0.25mm with roll groove;
11 03 0247	PSL404219	11mm Alu. Crimp Cap, clear lacquered, centre hole;	UltraClean Silicone white/PTFE red, 45° shore A, 1.3mm
11 03 0196	PSL404231	11mm Alu. Crimp Cap, clear lacquered, centre hole;	PTFE red/Silicone white/PTFE red, 45° shore A, 1.0mm
11 03 0464		11mm Alu. Crimp Cap, clear lacquered, centre hole;	Silicone white/PTFE blue, cross-slitted, 55° shore A, 1.5mm

Further 11mm Aluminium Crimp Seals with different Septa or with Crimp Caps in different colours are available!
Packaging Unit: 100 pieces per PE-Bag

Seals for Snap Ring ND11 (with assembled Septum)

11 15 1151	186000303 (alt.)	11mm PE Snap Ring Cap, blue, centre hole;	UltraClean Silicone white/PTFE red, 45° shore A, 1.3mm
11 15 1325	186002650 (alt.)	11mm PE Snap Ring Cap, red, centre hole;	UltraClean Silicone white/PTFE red, 45° shore A, 1.3mm
11 15 1269	186000304 (alt.)	11mm PE Snap Ring Cap, blue, centre hole;	Silicone white/PTFE blue, cross-slitted, 55° shore A, 1.0mm
11 15 1326	186002647 (alt.)	11mm PE Snap Ring Cap, red, centre hole;	Silicone white/PTFE blue, cross-slitted, 55° shore A, 1.0mm

Further 11mm Snap Ring Seals with different Septa resp. with red, transparent, green or yellow caps are available!
Packaging Unit: 100 pieces per PE-Bag

HPLC

4



Consumables for Waters

Waters Alliance®
Alliance® 2690/2695
Alliance® 2790/2795
ACQUITY UPLC™
Waters Alliance® HT System



We can also supply barcode labelled or pre-sealed vials!

Seals for Screw Neck ND10 (with assembled Septum)

10 15 1256	10mm PP Screw Cap, black, centre hole;	Natural Rubber red-orange/TEF transparent, 60° shore A, 1.3mm
10 15 1257	10mm PP Screw Cap, black, centre hole;	UltraClean Silicone white/PTFE red, 45° shore A, 1.3mm
10 15 0744	10mm PP Screw Cap, black, centre hole;	Silicone white/PTFE beige, 45° shore A, 1.5mm
10 15 1328	10mm PP Screw Cap, black, centre hole;	Silicone white/PTFE blue, with slit, 55° shore A, 1.5mm
10 15 1258	10mm PP Screw Cap, black, centre hole;	PTFE red/Silicone white/PTFE red, 45° shore A, 1.0mm
10 08 0742	10mm PP Screw Cap, black, centre hole;	
09 02 0723	9mm Septa: Silicone white/PTFE beige, 45° shore A, 1.5mm	
Packaging Unit: 100 pieces per PE-Bag (Seals + Screw Caps); 1000 pieces per PE-Bag (Septa)		

Shell Vials for WATERS WISP

Shell Vials, 1ml 96-Position Carousel			Shell Vials, 4ml 48-Position Carousel			
1:1 Drawings						
LPP Art.-No.	08 14 0513*	08 14 0595*	05 09 0706	15 14 0548	15 14 0562	
See Waters	WAT025054C	WAT025053C	alt. WAT072294	WAT025051	WAT025050	
Art.-No.						
LPP Art.-Description	1.0ml Shell Vial, 40 x 8.2mm, clear glass, 8mm PE-Plug, transparent, insertion barrier	1.0ml Shell Vial, 40 x 8.2mm, amber glass, 8mm PE-Plug, transparent, insertion barrier	0.1ml Micro-Insert, 34 x 5mm, clear glass, 13mm top	4.0ml Shell Vial, 44.5 x 14.65mm, clear glass, 15mm PE-Plug, transparent	4.0ml Shell Vial, 44.5 x 14.65mm, amber glass, 15mm PE-Plug, transparent	0.3ml Micro-Insert, 43.45 x 6mm, clear glass, 13mm top
LPP Packaging Unit	100 pieces per PP-Box (Vials) 100 pieces per PE-Bag (Plugs)	10 x 100 pieces per PP-Box	10 x 100 pieces per PP-Box	100 pieces per PP-Box (Vials) 100 pieces per PE-Bag (Plugs)	10 x 100 pieces per PP-Box	

*: alternatively available with a PE-Plug without insertion barrier for the Micro-Insert (==> 08 14 1168 resp. 08 14 1169). Thus the plug is not as sturdy and penetration and valve effect are further improved.

HPLC

5



Consumables for Waters

Waters WISP
GPC 2000
Waters 717



Screw Neck ND13 48-Position Carousel					
1:1 Drawings					
LPP Art.-No.	13 09 0222	13 09 1335	13 09 0280	13 09 1336	40 09 0146
See Waters	186000840	186001135	186001135	186001135	WAT015199
Art.-No.	PSL407498	PSL407498	PSL407499	PSL407499	WAT072708
LPP Art.-Description	4ml Screw Neck Vial, 45 x 14.7mm, clear glass, 1" hydrol. class	4ml Screw Neck Vial, 45 x 14.7mm, clear glass, 1" hydrol. Class, label + filling lines	4ml Screw Neck Vial, 45 x 14.7mm, amber glass, 1" hydrol. class	4ml Screw Neck Vial, 45 x 14.7mm, amber glass, 1" hydrol. Class, label + filling lines	0.3ml Micro-Insert, 40 x 6mm, clear glass, 15mm top
LPP Packaging Unit	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	100 pieces per PP-Box	10 x 100 pieces per PP-Box
					For Micro-Insert 40 09 0146

Seals for Screw Neck ND13 (with assembled Septum)

13 15 0456	PSL410208	13mm PP Screw Cap, black, centre hole;	Natural Rubber red-orange/TEF transp., 60° shore A, 1.3mm
13 15 0439	PSL410408	13mm PP Screw Cap, black, closed top;	Natural Rubber red-orange/TEF transp., 60° shore A, 1.3mm
13 15 0815	alt. 186000841	13mm PP Screw Cap, black, centre hole;	Silicone cream/PTFE red, 55° shore A, 1.5mm
13 15 0292		13mm PP Screw Cap, black, centre hole;	PTFE red/Silicone white/PTFE red, 45° shore A, 1.0mm
13 15 1293	alt. 186000842	13mm PP Screw Cap, black, centre hole;	Silicone white/PTFE blue, cross-slitted, 55° shore A, 1.5mm
13 08 0166	600000162	13mm PP Screw Cap, black, centre hole;	
12 02 0168	alt. WAT072711 (144 pcs.) WAT073005 (1440 pcs.) WAT072714 (144 pcs.)	12mm Septa: PTFE virginal, 53° shore D, 0.25mm	
		Further 13mm Srew Seals with different Septa resp. with closed top caps are available! Packaging Unit: 100 pieces per PE-Bag (Seals) 100 pieces per PE-Bag (Caps); 1,000 pieces per PE-Bag (Lines)	

HPLC

6



Consumables for Waters



Special La-Pha-Pack® 2in1 KITS for WATERS

SHORT THREAD ND9			
Article No.		Vial	+ Seal
11 24 1143	alt. 186000272C	11 09 0519 + 09 15 0838	
11 24 1628	alt. 186000272C	11 09 0500 + 09 04 1533	
11 24 1238	alt. 186000307C	11 09 0519 + 09 15 0869	
11 24 1226	alt. 186000307C	11 09 0500 + 09 04 1534	
11 24 1342	alt. 186000846C	11 09 0520 + 09 15 0838	
11 24 1...	186000846C	11 09 0520 + 09 04 1533	
11 24 1573	alt. 186000847C	11 09 0520 + 09 15 0869	
11 24 1...	186000847C	11 09 0520 + 09 04 1534	
11 24 1...	alt. 186000326C	11 09 0620 + 09 15 0838	
11 24 1...	186000326C	11 09 0620 + 09 04 1533	
11 24 1129	alt. 186000327C	11 09 0620 + 09 15 0869	
11 24 1...	186000327C	11 09 0620 + 09 04 1534	
11 24 1...	186001126C	11 09 0999 + 09 04 1533	
11 24 1...	186001128C	11 09 0999 + 09 04 1534	

4ml SHELL VIALS			
Article No.		Vial	+ Seal
15 27 1065	See. WAT025051	15 09 0539 + 15 08 0547	
15 27 1066	See WAT025050	15 09 0540 + 15 08 0547	

Further 2in1 Kits available upon request!

Others

08 08 0027	WAT210684	8mm PP Screw Cap, black, centre hole
08 15 0293	WAT094174	8mm PP Screw Cap, black, centre hole; Silicone white/PTFE red, 45° shore A, 1.3mm
08 02 0039	WAT058876	8mm Septa: Silicone white/PTFE blue, with slit, 55° shore A, 0.9mm
08 02 0103	WAT210685	8mm Septa: Silicone white/PTFE red, 45° shore A, 1.3mm
08 02 0177	WAT058886	8mm Septa: PTFE virginal, 53° shore D, 0.25mm
20 09 1405	WAT223176	10ml Headspace-Vial, 46 x 22.5mm, clear glass, 1 st hydrol. class, DIN-crimp neck, rounded bottom
20 03 0142	600000138	20mm Alu. Crimp Cap, plain, centre hole; Silicone transp. blue/PTFE white, 45° shore A, 3.0mm

SCREW NECK ND10			
Article No.		Vial	+ Seal
10 32 1612	WAT270946C	10 09 0743 + 10 15 0744	
10 32 1283		10 09 1197 + 10 15 1257	
10 32 1329		10 09 0743 + 10 15 1328	
10 32 1330		10 09 1196 + 10 15 1328	
10 32 1331		10 09 1197 + 10 15 1328	

PP MICRO-VIALS			
Article No.		Vial	+ Seal
11 24 1130	alt. 186002639	11 19 0932 + 09 15 0869	
11 24 1...	186002639	11 19 0932 + 09 04 1534	
11 24 1476	alt. 186002640	11 19 0932 + 09 15 0838	
11 24 1...	186002640	11 19 0932 + 09 04 1533	
11 24 1...	alt. 186002635	11 19 1205 + 09 15 0838	
11 24 1557	alt. 186002636	11 19 1205 + 09 15 0869	

SCREW NECK ND13			
Article No.		Vial	+ Seal + Septa
13 28 1069	alt. 186000838C	13 09 0222 + 13 15 0815	
13 28 1071	WAT073018	13 09 0222 + 13 08 0166 + 12 02 0168	
13 28 1541	alt. 186000839C	13 09 0222 + 13 15 1293	
13 28 1074	alt. 186001133C	13 09 0280 + 13 15 0815	
13 28 1076		13 09 0280 + 13 08 0166 + 12 02 0168	
13 28 1...	alt. 186001134C	13 09 0280 + 13 15 1293	

HPLC

Numerical Index / Cross References Waters - La-Pha-Pack

Waters Art.-No.	LPP Article-No.	These articles are included in the Waters brochure!
186000272C	11 24 1628 (alt)	
186000273	11 24 1143 (alt)	
186000273DV	11 09 0519	
186000274	09 04 1533	
186000303	09 04 1220 (alt)	
186000304	09 15 0838 (alt)	
186000305	11 15 1151 (alt)	
186000305	11 15 1289 (alt)	
186000305	09 04 1534	
186000326C	11 24 1... (11 09 0620 + 09 04 1533)	
186000326C	11 24 1... (alt.) (11 09 0620 + 09 15 0838)	
186000327C	11 24 1... (11 09 0620 + 09 04 1534)	
186000385	11 24 1129 (alt)	
186000386	13 28 1069 (alt)	
18600039C	13 28 1541 (alt)	
186000840	13 09 0222	
186000841	13 15 0875 (alt)	
186000842	13 15 1293 (alt)	
186000843	11 24 1... (11 09 0519 + 09 04 1533)	
186000844	11 24 1143 (alt)	
186000845	11 09 0519	
186000846C	11 24 1... (11 09 0520 + 09 04 1533)	
186000847C	11 24 1... (11 09 0520 + 09 04 1534)	
186000848	11 24 1573 (alt)	
186000848	11 09 0520	
186000848DV	11 09 1242	
186001126C	11 24 1... (11 09 0999 + 09 04 1533)	
186001126C	11 24 1... (11 09 0999 + 09 15 0838, alt)	
186001128C	11 24 1... (11 09 0999 + 09 04 1534)	
186001128C	11 24 1... (11 09 0999 + 09 15 0869, alt)	
186001133C	13 28 1074 (alt)	
186001134C	13 28 1... (alt.) (13 09 0280 + 13 15 1293)	
186001135	13 09 0280	
186002626	11 19 0932	
186002628	11 19 0933	
186002634	11 19 1205 (alt)	
186002635	11 24 1... (alt.) (11 19 1205 + 09 04 1533)	
186002635	11 24 1... (alt.) (11 19 1205 + 09 15 0838)	
186002636	11 24 1... (11 19 1205 + 09 04 1534)	
186002636	11 24 1557 (alt)	
186002639	11 24 1... (11 19 0932 + 09 04 1534)	
186002640	11 24 1643	
186002640	11 24 1476 (alt)	
186002642	11 26 1... (alt.) (11 19 0933 + 11 15 1151)	
186002641	11 26 1... (alt.) (11 19 0933 + 11 15 1269)	
186002647	11 15 1326 (alt)	
186002650	11 15 1325 (alt)	
186002802	11 09 0620	

Fields marked in yellow may include Art.-Nos. of 2in1 KITS that have not yet been issued in our computer system (e.g. 11 24 1...). However, we indicated instead Art.-No. of Vial + Seal.

Bold listed LPP Article Nos. are the most identical ones besides some other alternatives which, however, would suit from their functionality as well.

Other WATERS Cross references	
186002127	09 15 1...
186002128	09 15 1179 (alt)
186002129	09 15 1178 (alt)
186002130	09 15 1332 (alt)
PSL400123	20 09 0289
PSL403803	08 09 0606
PSL404219	08 03 0249
PSL407167	11 09 0415
PSL407498	13 09 0222
PSL407499	13 09 0280
PSL407500	13 08 0336
PSL408004	09 09 0426
PSL408008	08 03 0451

Waters Art.-No.	LPP Article-No.	These articles are included in the Waters brochure!
186002904	11 09 0999	
600000138	20 03 0142	
600000162	13 08 0166	
PSL404219	11 03 0247	
PSL404231	11 03 0196	
PSL408027	11 03 0339	
PSL410208	13 15 0456	
PSL410408	13 15 0439	
PSL412115	11 03 0300	
PSL904301	11 06 0006	
WAT015199	40 09 0146	
WAT0210685	08 02 0009 (alt.)	
WAT025050	15 14 0562	
WAT025050	15 27 1066 (see)	
WAT025051	15 14 0548	
WAT025051	15 27 1065 (see)	
WAT025053C	08 14 0595	
WAT025054C	08 14 0513	
WAT058874	10 15 0744	
WAT058875	10 08 0742	
WAT058876	08 02 0039	
WAT058886	08 02 0177	
WAT063300	10 09 0743	
WAT072294	05 09 0706 (alt.)	
WAT072706	50 13 0147	
WAT072710	13 09 0222	
WAT072711 (144)	13 08 0166 (alt.)	
WAT072714 (144)	12 02 0168	
WAT073005 (144)	13 28 1071	
WAT073916	08 09 0865	
WAT094170	06 09 1343	
WAT094174	08 15 0293	
WAT094174	08 15 0427 (alt.)	
WAT094219	11 09 0644	
WAT094220	11 09 0645	
WAT094222	11 09 0356	
WAT094223	11 09 0477	
WAT210684	08 08 0027	
WAT210685	08 02 0103	
WAT223176	20 09 1405	
WAT270946C	10 32 1612	

Other WATERS Cross references	
PSL408026	11 02 0176
PSL408027	11 01 0173
PSL408028	08 09 0305
PSL408067	20 09 0297
PSL408068	20 09 0342
PSL408101	20 01 0408
PSL408153	13 15 0862
PSL408225	08 09 0276
PSL408234	20 03 0030
PSL408789	11 19 1022
PSL412109	08 03 0268
PSL414168	20 02 0206
PSL904303	20 06 0008
WAT078515	08 08 0488

All data indicated herein has been made at best knowledge. We do not take over any guarantees for correctness nor for completeness. Status 2005