


Yarra™ 5 µm SEC/GFC PREP

new Higher Performance for Preparative BioSeparations at a Lower Price

Enjoy the same selectivity and ultra-high efficiency of Yarra 3 µm for your preparative gel filtration applications. Yarra SEC PREP features a 5 µm particle size version of the original Yarra 3 µm particle with the same chemistry on a 21.2 mm ID column for preparative purification, desalting, and characterization of biomolecules. Yarra 5 µm PREP is available at a price and performance that other preparative silica GFC columns cannot match.

Yarra 3 µm SEC-2000 vs. Yarra 5 µm SEC-2000 PREP Column



Yarra 5 µm PREP SEC Columns (mm)	Preparative	SecurityGuard™ Cartridges (mm)
Phases	300 x 21.2	15 x 21.2**
		/ea
Yarra 5 µm SEC-2000 PREP	00H-4619-P0	AJ0-8588
Yarra 5 µm SEC-3000 PREP	00H-4620-P0	AJ0-8589
Yarra 5 µm SEC-4000 PREP	00H-4621-P0	AJ0-8950

**PREP SecurityGuard™ Cartridges require holder, Part No.: AJ0-8223 for ID: 18 - 29 mm